

Prospect Designs Supine Stander



The Prospect Designs **Supine Stander** is a versatile, easy to adjust option for children who need supported upright weight bearing and for whom prone standing is not appropriate. The child is always positioned at the lowest possible height regardless of their size and the stander “grows” up with the child. The wide range of angle adjustment allows the stander to be positioned nearly horizontal while placing the child in the stander. The stander can then be easily raised to the desired angle for the child to bear weight.



Available in two sizes: (SS1) (SS2)

FEATURES

-Basic stander includes: two pairs of lateral supports (for the trunk and the hips), one 2" velcro strap for either hip or trunk laterals, sectional back pads, posterior knee pads with padded strap, shoe holders and locking casters

SPECIFICATIONS

-Footrest to top of back board	27" to 36	35" to 45"
-Range of recline angle	90° to 5°	90° to 25°
-Width between laterals	6 1/2" to 11"	8" to 13"
-Width of main supine board	15"	17"
-Base--width X length (fully extended)	24" X 36	27" X 36"



Small Supine Stander (SS1)
with Contoured Headrest, Chest
Harness and clip on tray



Large Supine Stander (fully reclined) (SS2) with
Contoured Headrest, Chest Harness, and Lateral
Knee supports

ACCESSORIES



Lateral kneesupports
Adjustable for height and width. Can be used inside or outside the knee
LKS

Shoulder Protractors

45° support pads mounted in the same slots as the standard lateral pads help to block unwanted upper extremity extension or to keep arms more mid-line for children with low tone. SHPRS (Small) is used on the SS1 SHPRL (large) is used on the SS2



SHPRS, SHPRL

Adjustable armrests

Adjustable armrests can be added to the basic supine stander. They are angle adjustable and may be flipped back completely out of the way. They are height adjustable by removing bolts and repositioning them in a series of holes on the sides of the chair.



AR

Slide On Tray

The tray slides onto the adjustable armrests and is secured with two spring pins that engage in a series of holes on the armrests.



TRS

Clip on tray

Attaches to the chair with two spring clips that engage in a series of holes on the sides of the chair. The height is determined by the choice of holes. The angle is adjusted by supports attached to the tray.



TRC



Chest harness

Four point bib style harness. Secured with 4 fastex buckles on 1" webbing. Two sizes correspond to the two sizes of stander.

CHS1, CHS2

Headrest and Headrest Mounting Options



Contoured headrest pad with adjustable height and depth
CH-ADH

Flat headrest pad and laterals with adjustable height
FHL-AH



Contoured headrest pad with adjustable height
CH-AH

Flat headrest pad and laterals with adjustable height and depth
FHL-ADH



Flat headrest pad with adjustable height and depth
FH-ADH

Flat headrest pad with adjustable height
FH-AH

