



- Lightweight, easy-to-handle and set up
- Folds in half, carries like a suitcase*
- Durable welded construction
- Full platform provides excellent stability
- Anti-slip, high traction surface
- Accommodates wheelchairs and scooters with various wheel configurations
- Safety DVD and steel security pins included



Custom widths are available. Give your specifications to your local PVI retailer!

Item Number	Length	Width	Folded Size	Weight	Capacity**
SFW230*	2'	30"	24" x 16" x 4 $\frac{7}{8}$ "	10 lbs.	800 lbs.
SFW330	3'	30"	36" x 16" x 4 $\frac{7}{8}$ "	15 lbs.	800 lbs.
SFW430	4'	30"	48" x 16" x 4 $\frac{7}{8}$ "	20 lbs.	800 lbs.
SFW530	5'	30"	60" x 16" x 4 $\frac{7}{8}$ "	24 lbs.	800 lbs.
SFW630	6'	30"	72" x 16" x 4 $\frac{7}{8}$ "	28 lbs.	800 lbs.
SFW730	7'	30"	84" x 16" x 4 $\frac{7}{8}$ "	34 lbs.	800 lbs.
SFW830	8'	30"	96" x 16" x 4 $\frac{7}{8}$ "	46 lbs.	800 lbs.

Not sure which ramp length you need? Refer to pages 18 - 19 to help you.



*Handle not included on SFW230. ** Maximum ramp capacity: 400 lbs. (one axle) – 800 lbs. (two axles).