



Stores Upright in Garage



Simply Wheel Out to Required Location

- Closure strap locks ramp panels together
- Easy-to-handle & set up
- Durable welded construction
- Full platform provides excellent stability
- Anti-slip, high traction surface
- Accommodates wheelchairs and scooters with various wheel configurations
- Stores upright in garage
- Safety DVD and steel security pins included

### Effortless Wheel Design



Quick Release Pins used to secure Latch Straps



WAB1030 &amp; WAB1230 Feature 2 3/4" Long Top Hook



Solid / Stable Surface

LIMITED  
**LIFETIME  
WARRANTY**

| Item Number | Top Hook Length | Length | Width | Folded Size             | Weight  | Capacity* |
|-------------|-----------------|--------|-------|-------------------------|---------|-----------|
| WAB1030     | 2 3/4"          | 10'    | 30"   | 66 1/4" x 36 3/4" x 16" | 84 lbs. | 600 lbs.  |
| WAB1230     | 2 3/4"          | 12'    | 30"   | 78 1/4" x 36 3/4" x 16" | 98 lbs. | 600 lbs.  |
| WABU1030    | 9 3/4"          | 10'    | 30"   | 63 1/4" x 36 3/4" x 16" | 83 lbs. | 600 lbs.  |
| WABU1230    | 9 3/4"          | 12'    | 30"   | 75 1/4" x 36 3/4" x 16" | 97 lbs. | 600 lbs.  |

Not sure which ramp length you need? Refer to pages 2 - 3 to help you.

\* Maximum ramp capacity: 300 lbs. (one axle) – 600 lbs. (two axles).



**PVI**