

  
**QUICKIE®**

**Xtender™**

Allows you to  
go the distance!



## Live Without Limits™

The Quickie® Xtender™ provides power assist to a manual wheelchair, allowing users to maintain the look, control, and functionality of a manual wheelchair. The Xtender increases the distance traveled and the efficiency of propulsion without additional physical exertion.



da Vinci  
AWARDS®

### Powered by Yamaha

- Motorized, power assist wheels increase the force applied to the handrims by four times
- Electronically synchronized units compensate for arm strength imbalance for optimal response and straight propulsion
- Advanced lightweight battery technology includes a microprocessor for increased performance
- Features downhill speed control and rollback assistance for improved performance during ascents and descents
- Available with nickel or lithium battery options

### Increased Access and Independence

- Reduces energy expenditure by as much as 45%, reduces user fatigue, and increases range of mobility
- Provides excellent handrim grip for users with limited hand function
- May decrease risk of long term shoulder and wrist injuries

### Ease of Use

- Quick-release wheels allow for easy transitions between manual and assisted propulsion
- No programming required
- Easy to transport and store
- Bright LED light alerts user when unit is on
- Battery level lights indicate remaining battery power

SPECIFICATIONS	Quickie Xtender
HCPCS Code	E0986
Available On	On Chair: Quickie 2, 2HP, Q7, GT, GP Series Through Built-4-Me: Quickie 2 Lite, QRi
Torque Multiplier	1:4
Overall Weight	with Lithium Battery: 36 lbs. off chair with Nickel Battery: 34 lbs. off chair
Weight Capacity	265 lbs.
Battery Options	Lithium or Nickel
Battery Capacity	Lithium: 9 hours Nickel: 7 hours
Time to Full Battery Charge	Lithium: 3 to 4.5 hours Nickel: 2.5 to 3 hours
Average Range	Lithium: 17 miles Nickel: 9.5 miles
Available Tires	Pneumatic
Available Handrims	Stainless Steel or Plastic Coated