



# Jaz

## Mobile Shower Commode Chairs



Jaz-SP Self Propel  
(shown with optional  
SoftJaz seat overlay)

Jaz mobile shower commode chairs from Raz Design offer exceptional comfort and function at an affordable price ... but with standard features that are normally found on more expensive chairs!



Jaz-AP  
Attendant Propel

Jaz mobile shower commode chairs are exceptionally durable and easy to clean, ideal for infection-control environments.


FEATURES



Molded, composite frame – durable and easy to clean



Height-adjustable without tools – for toilet bowl clearance and matching transfer heights




Contoured seat for best-in-class comfort



Swing-away, removable foot supports – facilitates transfers



Soft and flexible back support – comfortable and easy to clean



Flip-up arm supports – aids in transfers



Dual-locking casters (swivel and roll) – stabilize chair for toileting, bathing and transfers



Commode pan with lid – includes anti-splash baffle and high-gloss surface to ease cleaning

OPTIONS



SoftJaz seat overlay – molded and contoured for pressure distribution and unrivaled comfort



Pelvic positioning belt with side-release buckle



Chest belt – 1-piece with hook and loop fastener



Snap-on Bodypoint Calf strap – dries quickly



Arm support locks – provide more security during transfers



Swing-away laterals – provide lateral trunk support



Head support - molded pad with fixture, adjusts in height and fore/aft



Quick-release axles (Jaz-SP) – makes transporting the chair easier

SPECS

SPECIFICATIONS	Jaz-AP	Jaz-SP
WEIGHT CAPACITY	150 kg (330 lbs)	150 kg (330 lbs)
SEAT WIDTH	48 cm (19")	48 cm (19")
SEAT DEPTH (front of back post to front of seat)	43 cm (17")	43 cm (17")
BACK WIDTH	46 cm (18")	46 cm (18")
BACK HEIGHT (from top of seat to top of upholstery)	43 cm (17")	43 cm (17")
TOILET HEIGHT CLEARANCE:	39-51 cm (15.4"-20.4")	39 cm-51 cm (15.4"-20.4")
TOILET WIDTH CLEARANCE	40 cm (15.7")	39 cm (15.4")
ARM SUPPORT HEIGHT ABOVE SEAT*	20 cm (8")	20 cm (8")
SEAT-TO-FLOOR HEIGHT* – 2.5 cm (1") increments	48-61 cm (19"-24")	48-61 cm (19"-24")
OVERALL WIDTH	58 cm (23")	66 cm (26")
PRODUCT WEIGHT	17 kg (38 lbs)	22 kg (48 lbs)
FRAME MATERIAL	Molded Composite	Molded Composite

\* without SoftJaz seat overlay