



Raz-SP

MSRP – USA – USD
Effective July 1, 2020

Date: _____

Dealer / Location: _____ T: _____ F: _____

PO# / Client: _____ Contact: _____ Email: _____

Raz-SP Self Propel Mobile Shower Commode Chair. Includes height-adjustable stainless steel frame, flip-up padded arm supports, commode pan, foot supports, 22" or 24" rear wheels and two 5" dual-locking casters. Standard seat-to-floor height range is 20"-24" Exclusive to Raz, all seats include the Ischial Pelvic Alignment System [IPAS], which allows for 2" fore/aft adjustment. **User weight capacity 400 lbs.** Please enter special set-up requests, such as arm support height, in Notes [Sec.16].

For custom modifications of frames and accessories, please contact Raz Design. To assist in evaluation and justification see our [Evaluation and Justification Form](#).

1. Chair Codes by size – Must Choose One

<input type="checkbox"/> Z20018	Raz-SP w/ 18"W frame – Standard frame size, will fit over most toilets ¹	\$1,950
<input type="checkbox"/> Z200__*	Non-Standard Width Raz-SP – Must add upcharge for custom frame width: see below * add custom width to complete code (16/20/22) – Adjusta-Back must be selected for all custom frame widths [Sec.5]	\$1,950
<input type="checkbox"/> Z2SW16	Upcharge for custom 16"W frame ² – will not fit over most toilets	\$495
<input type="checkbox"/> Z2SW20	Upcharge for custom 20"W frame ³	\$495
<input type="checkbox"/> Z2SW22	Upcharge for custom 22"W frame ⁴	\$595

¹ Fits all 18"W seats

² Visco-Foam Interface seats can be modified to fit.

³ 18"W seats can be selected but must also select Seat Support Adapter Kit Z4SSKAS [Sec.2]

⁴ Z4V23N Visco Foam Interface HD seat can be modified to fit. Seat depth range (IPAS) will be 18"-20".

Custom Base Frame Length and Height

<input type="checkbox"/> Z2LNGF3	Upcharge – custom 3" Longer Frame recommended for seat depths 20"/21", req for seat depths 21"-24"	\$395
<input type="checkbox"/> Z2LOWF4	Upcharge for custom 4" Lower Frame – reduces seat height range to 16"-18" from standard height range of 20"-24" with molded seat – To lower chair by an additional 1" (15"-17"), add 4" caster upgrade option Z64SP [Sec.7] – Compatible only with 22" wheels [Sec.8]	\$495

2. Seats – Standard

Molded Seats (thickness ~2.5"-3")			Visco-Foam Interface Seats (thickness ~3.5")		
<input type="checkbox"/> Z405	18"W x 16"D – Front Access	Standard	<input type="checkbox"/> Z436BRI	18"W x 16"D – Front Bridged	\$425
<input type="checkbox"/> Z410	18"W x 16"D – Front Bridged	N/C	<input type="checkbox"/> Z436FRO	18"W x 16"D – Front Access	\$425
<input type="checkbox"/> Z415	18"W x 19"D – Front Access	\$135	<input type="checkbox"/> Z436LTO	18"W x 16"D – Left Access	\$425
<input type="checkbox"/> Z420	18"W x 19"D – Front Bridged	\$135	<input type="checkbox"/> Z436RTO	18"W x 16"D – Right Access	\$425
Specialty Seats			<input type="checkbox"/> Z436RRO	18"W x 16"D – Rear Access	\$425
<input type="checkbox"/> Z4SOM	Shower-Only Seat – 18"W x 16"-21"D	\$395	<input type="checkbox"/> Z439BRI	18"W x 19"D – Front Bridged	\$495
<input type="checkbox"/> Z4EJRM	E&J Replica Seat – Molded – 18"W x 17"D	\$295	<input type="checkbox"/> Z439FRO	18"W x 19"D – Front Access	\$495
<input type="checkbox"/> Z4EJRV	E&J Replica Seat – Visco Foam – 18"W x 17"D	\$895	<input type="checkbox"/> Z439LTO	18"W x 19"D – Left Access	\$495
<input type="checkbox"/> Z4V23N	23"W x 18"D HD – Front Access Notched at rear to maintain fore/aft adjustment range (IPAS) – 22"W frame	\$695	<input type="checkbox"/> Z439RTO	18"W x 19"D – Right Access	\$495
			<input type="checkbox"/> Z439RRO	18"W x 19"D – Rear Access	\$495
			<input type="checkbox"/> Z4SSKAS	Seat Support Kit for Raz-AP/SP – for seat 2" narrower than seat frame	\$145
			<input type="checkbox"/> Z4GSB__*	Custom Gimballed Seat Brackets – for adding dump angle (pr) Specify additional height at front of seat to complete code (1.0"/1.5"/2.0")	\$295

3. Seats – Custom

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required.

<input type="checkbox"/> Z4CS	Custom Seat (enter specs below) Select if both width and depth are ≤ 21"	\$795
<input type="checkbox"/> Z4CEJR	Custom E&J Replica Seat (enter specs below) Select if both width and depth are ≤ 21"	\$995
<input type="checkbox"/> Z4CSHD	Heavy Duty Custom Seats (enter specs below) Select if width and/or depth is ≥ 22"	\$1,095

Custom Specifications (for above custom seats)

Seat Frame Width:	-- 16" / 18" / 20" / 22"
Custom Seat Width:	-- 16" / 18" / 20" / 22" / 24" (cannot be more than 2" wider than seat width)
Custom Seat Depth:	-- 15" / 16" / 17" / 18" / 19" / 20" / 21" / 22" / 23" / 24" [≥ 21" requires Long Frame (Z2LNGF3)]
Access Opening:	-- FR – Front / LT – Left / RT – Right / RR – Rear / BR – Bridged
Access Opening Width:	-- 2.5" / 3.0" / 3.5" / 4.0" / 4.5" / 5.0" (3.5"W is standard)
Aperture Width:	-- 3.0"-8.0" (6.5" recommended) – enter "NA" if aperture is not required

- Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seats.
- Custom Seat width should not be narrower than frame width.
- Custom Seat can be up to 2" wider than seat frame width to provide additional thigh support without ordering wider chair frame.
- The Ischial Pelvic Alignment System (IPAS) allows 2" fore/aft adjustment of the seat (seat depth and aperture position) on all seats.
- Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, etc.

4. Fixed Back Frames (18"W only)

These are called "Fixed" Backs because the arm supports are not adjustable in height. Includes tension-adjustable back upholstery, fixed-height flip-up padded arm supports, back angle adjustment (85°/90°/95°) with cane-style posts. Back heights are nominally measured from top of standard molded seat to top of back upholstery. Back upholstery can be positioned 2.5" higher than nominal height indicated below. e.g. 20"H can be increased to 22.5"H

<input type="checkbox"/> Z505	Fixed Back – 17" Height ¹	N/C	<input type="checkbox"/> Z506	Fixed Back – 20" Height ¹	\$70
-------------------------------	--------------------------------------	-----	-------------------------------	--------------------------------------	------

¹ Armrest Lock Kit (Z710) [Sec.7] is not included.

5. Adjusta-Back Frames

Includes height- and angle-adjustable arms with Armrest Lock Kit (Z710), tension-adjustable back upholstery, flip-up padded arm supports, back angle adjustment (85°/90°/95°). Back heights of 15"/17"/20" include cane-style posts, on which back upholstery can be positioned 2.5" higher than nominal height indicated below. e.g. 20"H can be moved up to 22.5"H. The 24"H Adjusta-Back includes straight posts, on which the upholstery cannot be raised.

Standard Width Adjust-Backs (18"W)			Non-Standard Width Adjust-Backs		
<input type="checkbox"/> Z50918	Adjusta-Back – 15"H	\$225	* Complete Adjusta-Back order code by adding desired back frame width (16"/20"/22") Note: 18" is the standard width.		
<input type="checkbox"/> Z51018	Adjusta-Back – 17"H	\$225	<input type="checkbox"/> Z509_*	Adjusta-Back – 15"H non-std width	\$225
<input type="checkbox"/> Z51518	Adjusta-Back – 20"H	\$295	<input type="checkbox"/> Z510_*	Adjusta-Back – 17"H non-std width	\$225
<input type="checkbox"/> Z52018	Adjusta-Back – 24"H	\$345	<input type="checkbox"/> Z515_*	Adjusta-Back – 20"H non-std width	\$295
<input type="checkbox"/> Z5DTL	Dual Triggers – both on left side ¹	N/C	<input type="checkbox"/> Z520_*	Adjusta-Back – 24"H non-std width	\$345
<input type="checkbox"/> Z5DTR	Dual Triggers – both on right side ¹	N/C	<input type="checkbox"/> Z515SPR_*	Reclining Adjusta-Back 30° recline, 20" back height with 2-Point Reclining Arms (incl Rear Anti-tippers Z66SPR)	\$895
<input type="checkbox"/> Z515SPR18	Reclining Adjusta-Back 30° recline, 20" back height with 2-Point Reclining Arms (incl Rear Anti-tippers Z66SPR)	\$895	<input type="checkbox"/> Z520SPR_*	Reclining Adjusta-Back 30° recline, 24" back height with 2-Point Reclining Arms (incl Rear Anti-tippers Z66SPR)	\$995
<input type="checkbox"/> Z520SPR18	Reclining Adjusta-Back 30° recline, 24" back height with 2-Point Reclining Arms (incl Rear Anti-tippers Z66SPR)	\$995			

¹ Dual Tilt Triggers option positions both triggers on one side, allowing for single-handed activation. Available for all back widths.

6. Custom Back Frame Width Narrower Than Chair Frame Width – Must choose back frame [Sec.5] (configurations below are not available for recliner backs)

NOTE: For back frame width variances on recliner frames, please request custom quotation.

<input type="checkbox"/> Z5NB__ __* Custom Back Frame 2" Narrower than base frame * Enter base and back frame widths to complete code. eg. Z5NB2018 specifies a 20"W chair frame with an 18"W back frame	\$345	Base Frame Width	16" ¹	18" ²	20" ³	22" ⁴
		Back Frame Width	14"	16"	18"	20"

Consider also using Arm Spacer Kit (Z718/Z719) as another solution for providing more width between the arms.



¹ Requires Custom Seat

² Compatible with all 18"W seats but with reduced seat depth range. Also compatible with Custom or modified Visco-Foam Interface seats.

³ Custom Seat recommended. Also compatible with all 18"W seats with Seat Support Kit Z4SSKAS [Sec.2].

⁴ Custom Seat recommended. Also compatible with modified Z4V23N (18" seat depth with limited adjustability)

7. Back Support Options – These back supports replace the standard, tension-adjustable upholstery and are recommended for infection-control environments.

<input type="checkbox"/> Z5S14 14"W Symphony Back ¹ – recessed solid back \$395	<input type="checkbox"/> Z749 Harmony Back ² (contoured solid back – 18"W x 13"H) Fits 17" or 20" back heights, 18"W back widths only \$245
<input type="checkbox"/> Z5S16 16"W Symphony Back ¹ – recessed solid back \$395	
<input type="checkbox"/> Z5S18 18"W Symphony Back ¹ – recessed solid back \$395	
<input type="checkbox"/> Z5S20 20"W Symphony Back ¹ – recessed solid back \$495	¹ Not compatible w/ 15" back height – Symphony Back panels are 13"H
<input type="checkbox"/> Z5S22 22"W Symphony Back ¹ – recessed solid back \$595	² Not compatible w/ 15" and 24" back heights. Not recommended for reclining backs.

8. Rear Wheel and Caster Options

<input type="checkbox"/> Z612R 22" Diameter Wheels w/ Textured Handrims Standard	<input type="checkbox"/> Z660 Quick-Release Axles (pr) \$95
<input type="checkbox"/> Z614R 24" Diameter Wheels w/ Textured Handrims N/C	<input type="checkbox"/> Z663 Wheel Lock Extension Kit (pr) \$55
<input type="checkbox"/> Z622VHR Vinyl-Coated Handrims for 22" (pr) \$265	<input type="checkbox"/> Z674 User-Actuated Front Directional Caster Locks w/ extension handles (pr) \$345
<input type="checkbox"/> Z624VHR Vinyl-Coated Handrims for 24" (pr) \$265	<input type="checkbox"/> Z64SP 4" Caster Upgrade – affects STF range (2x) \$140
<input type="checkbox"/> Z622PHR Vertical Projection Handrims – 22" (pr) \$465	
<input type="checkbox"/> Z624PHR Vertical Projection Handrims – 24" (pr) \$495	

9. Arm Supports and Accessories

<input type="checkbox"/> Z5BAR13	Flip-Up Arm Supports – 13”L pad (pr)	Std	<input type="checkbox"/> Z7PAML	Pivoting Arm Mount ³ – Left (for Z7FL/Z7TR)	\$165
<input type="checkbox"/> Z5BAR10	Flip-Up Arm Supports – 10”L pad (pr) (not compatible with Z7FS/Z7FL/Z7TR / Z7PAM)	N/C	<input type="checkbox"/> Z7PAMR	Pivoting Arm Mount ³ – Right (for Z7FL/Z7TR)	\$165
<input type="checkbox"/> Z710	Arm Support Lock Kit (for Z505/Z506) (pr)	\$75	<input type="checkbox"/> Z718	Arm Support Spacer Kit – 1/2” per side (for Z509/Z510/Z515/Z520) (pr)	\$60
<input type="checkbox"/> Z52PAS	2-Point Flip-up Arm Supports ¹ (pr) (standard with reclining backs)	\$345	<input type="checkbox"/> Z719	Arm Support Spacer Kit – 1” per side (for Z509/Z510/Z515/Z520) (pr)	\$115
<input type="checkbox"/> Z7FLL	Long Molded Flat Armpad ² – Left (14”L x 3.25”W, fits Z5BAR13)	\$95	<input type="checkbox"/> Z7PHG	Pivoting Hand Grips (pr)	\$395
<input type="checkbox"/> Z7FLR	Long Molded Flat Armpad ² – Right (14”L x 3.25”W, fits Z5BAR13)	\$95	<input type="checkbox"/> Z7APS_ _*	Anterior Posture Support Bar * Enter back frame width to complete code (14”/16”/18”) ^{3, 4, 5}	\$195
<input type="checkbox"/> Z7FSL	Short Molded Flat Armpad – Left (10”L x 3”W, fits Z5BAR10)	\$95	<input type="checkbox"/> Z7APS_ _*	Anterior Posture Support Bar * Enter back frame width to complete code (20”/22”) ^{3, 4, 5}	\$225
<input type="checkbox"/> Z7FSR	Short Molded Flat Armpad – Right (10”L x 3”W, fits Z5BAR10)	\$95			
<input type="checkbox"/> Z7TRL	Molded Arm Trough ² – Left (18.5”L x 4.5”W, fits Z5BAR13)	\$110			
<input type="checkbox"/> Z7TRR	Molded Arm Trough ² – Right (18.5”L x 4.5”W, fits Z5BAR13)	\$110			

¹ For chairs with same back width and frame width only
² May require Z7PAM for arm support to flip up and/or for lateral supports to swing away
³ Not compatible with Z5BAR10
⁴ Not compatible with Z7FL/Z7FS/Z7TR
⁵ Special length required with Z718 (+1”) and Z719 (+2”)

10. Foot and Leg Supports – Standard Foot Supports (Z705N) are included with the chair. The selection of any Foot / Leg Support Option below automatically replaces the relevant Standard Foot Support. e.g. By selecting Elevating Leg Support – Left (Z7ELRL), the chair will be supplied with one Foot Support (Z705N) on the right side and a Z7ELRL on the left side.

<input type="checkbox"/> Z705N	Foot Supports 16”-20.5” (pr)	Std	<input type="checkbox"/> Z7ELRL	Elevating Leg Support ⁴ – Left (includes Z7ADM)	\$475
<input type="checkbox"/> Z705N16	Foot Supports – 16”W frames 16”-20.5” (pr)	N/C	<input type="checkbox"/> Z7ELRR	Elevating Leg Support ⁴ – Right (includes Z7ADM)	\$475
<input type="checkbox"/> Z7FET2	Footrest Extension Tubes 2” longer ¹ (pr)	\$50	<input type="checkbox"/> Z7ELROFL	Forward Offset Elevating Leg Support ⁵ – Left (includes Z7ADM)	\$525
<input type="checkbox"/> Z7MFX8	MFX8 Foot Supports 8”-12” (pr)	\$75	<input type="checkbox"/> Z7ELROFR	Forward Offset Elevating Leg Support ⁵ – Right (includes Z7ADM)	\$525
<input type="checkbox"/> Z7MFX12	MFX12 Foot Supports 12”-16.5” (pr)	\$175	<input type="checkbox"/> Z7RLSL	Universal Residual Limb Support – Left	\$395
<input type="checkbox"/> Z7ADM	Adj Angle/Depth Molded Footplates ² (pr)	\$190	<input type="checkbox"/> Z7RLSR	Universal Residual Limb Support – Right	\$395
<input type="checkbox"/> Z7FET2A	Footrest Extension Tubes for Z7ADM (pr) (2” longer)	\$50			
<input type="checkbox"/> Z7ADABS	Angle/Depth Adj ABS Footplates (pr) for 16”W frames	\$190			
<input type="checkbox"/> Z7VS	V-Style Foot Support ³ (1-piece)	\$195			
<input type="checkbox"/> Z7FBAS	Flip-Back Foot Support ³ (1-piece)	\$195			

¹ Not compatible with Z7MFX8/Z7ADM/Z7VS/Z7FBAS
² Not compatible with Z20016 or Z7MFX8
³ Custom required for sizes other than 18”W
⁴ Z7ELRL/R compatible with frame widths ≥ 18” and adjusted seat depths ≤ 18” or ≤ 21” w/ 3” longer seat frames (Z2LNGF3)
⁵ Z7ELROFL/R compatible with frame widths ≥ 18” for adjusted seat depths 19”-20” on standard frame or 22”-23” w/ 3” longer seat frames (Z2LNGF3)

11. Foot and Leg Support Accessories

<input type="checkbox"/> Z7FSLK	Foot Support Lock (pr)	\$95	<input type="checkbox"/> Z7NCS	Neoprene Covers for Standard Footplates (pr)	\$90
<input type="checkbox"/> Z7CS_*	Snap-On Calf Strap – Bodypoint Aeromesh * Enter frame width to complete code (16"/18"/20"/22")	\$95	<input type="checkbox"/> Z7NC16	Neoprene Covers for narrow molded footplates (pr)	\$90
<input type="checkbox"/> Z737CPL	Calf Panel – Bodypoint Aeromesh ¹ Fits frame widths (16"/18"/20")	\$175	<input type="checkbox"/> Z7NCA	Neoprene Covers for Adj Angle/Depth Molded Footplates (pr)	\$90
<input type="checkbox"/> Z706	H-Strap – for 18"/20"W frames – hook and loop attachment	\$90	<input type="checkbox"/> Z7NCABS	Neoprene Covers for ABS Adj Angle/Depth Footplates Z7ADABS (pr)	\$90
<input type="checkbox"/> Z712	H-Strap – for 22"W frames – hook and loop attachment	\$110	<input type="checkbox"/> Z7OFLM	Lateral Offset Foot Support Receiver – moves foot support 2" laterally or medially Quantity (1 or 2) __	\$160
<input type="checkbox"/> Z708ADJ	Adjustable Heel Loops (pr) – not available for Z7VS / Z7FBAS	\$60	¹ Not compatible with Z7MFX8 / Z7FBAS / Z7VS		

12. Positioning Pads

<input type="checkbox"/> Z7HSASS	Head Support w/ Small Pad (10" x 4")	\$495	<input type="checkbox"/> Z7LGAS1L	Leg Adductor w/ Small Pad ¹ – Left	\$320
<input type="checkbox"/> Z7HSASL	Head Support w/ Large Pad (12" x 6")	\$545	<input type="checkbox"/> Z7LGAS1R	Leg Adductor w/ Small Pad ¹ – Right	\$320
<input type="checkbox"/> Z7OHIP	Lateral Offset Head Support Interface Plate– up to 22° of adjustment	\$95	<input type="checkbox"/> Z7LGAS2L	Leg Adductor w/ Large Pad ¹ – Left	\$320
<input type="checkbox"/> Z7LTS1L	Lateral Thoracic Support w/ Small Pad ^{1,2} – Left	\$320	<input type="checkbox"/> Z7LGAS2R	Leg Adductor w/ Large Pad ¹ – Right	\$320
<input type="checkbox"/> Z7LTS1R	Lateral Thoracic Support w/ Small Pad ^{1,2} – Right	\$320	<input type="checkbox"/> Z7HP1L	Hip Adductor w/ Small Pad ^{1,3} – Left	\$320
<input type="checkbox"/> Z7LTS2L	Lateral Thoracic Support w/ Large Pad ^{1,2} – Left	\$320	<input type="checkbox"/> Z7HP1R	Hip Adductor w/ Small Pad ^{1,3} – Right	\$320
<input type="checkbox"/> Z7LTS2R	Lateral Thoracic Support w/ Large Pad ^{1,2} – Right	\$320	<input type="checkbox"/> Z7HP2L	Hip Adductor w/ Large Pad ^{1,3} – Left	\$320
<input type="checkbox"/> Z742L	Lateral/Adductor Extension Plate (5.5" x 3.5") – Left	\$25	<input type="checkbox"/> Z7HP2R	Hip Adductor w/ Large Pad ^{1,3} – Right	\$320
<input type="checkbox"/> Z742R	Lateral/Adductor Extension Plate (5.5" x 3.5") – Right	\$25	¹ Adjustable Pad Sizes: Small – 5.5"x 3.5" / Large – 6.75" x 4.3".		
			² Pads on Lateral Thoracic Supports can be adjusted medially up to 3" per side. For additional medial offset, specify Z742 extension plate.		
			³ Fits Adjusta-Back only		

13. Pelvic and Chest Belts – A chest belt should not be used without a pelvic belt

Pelvic Belts				Chest Belts					
<input type="checkbox"/>	ZPBS	Pelvic Belt – Small (frame widths 14"-16", 2-piece w/ side-release buckle)	\$35	<input type="checkbox"/>	ZCBS	Chest Belt ² – Small (frame widths 14"-16", 1-piece w/ hook and loop closure)	\$35		
<input type="checkbox"/>	ZPBM	Pelvic Belt – Medium (frame widths 16"-18", 2-piece w/ side-release buckle)	\$35	<input type="checkbox"/>	ZCBM	Chest Belt ² – Medium (frame widths 16"-18", 1-piece w/ hook and loop closure)	\$35		
<input type="checkbox"/>	ZPBL	Pelvic Belt – Large (frame widths 18"-22", 2-piece w/ side-release buckle)	\$45	<input type="checkbox"/>	ZCBL	Chest Belt ² – Large (frame widths 18"-22", 1-piece w/ hook and loop closure)	\$45		
<input type="checkbox"/>	ZPBXL	Pelvic Belt – XLarge (frame widths 22"-30", 2-piece w/ side-release buckle)	\$55	<input type="checkbox"/>	ZCBXL	Chest Belt ² – XLarge (frame widths 22"-30", 1-piece w/ hook and loop closure)	\$55		
<input type="checkbox"/>	ZPBICS	Infection Control Pelvic Belt ¹ – Small (frame widths 14"-16", 2-piece w/ side-release buckle)	\$75	<input type="checkbox"/>	ZCBICS	Infection Control Chest Belt ^{1,2} – Small (frame widths 14"-16", 1-piece w/ side-release buckle)	\$75		
<input type="checkbox"/>	ZPBICM	Infection Control Pelvic Belt ¹ – Medium (frame widths 16"-18", 2-piece w/ side-release buckle)	\$85	<input type="checkbox"/>	ZCBICM	Infection Control Chest Belt ^{1,2} – Medium (frame widths 16"-18", 1-piece w/ side-release buckle)	\$75		
<input type="checkbox"/>	ZPBICL	Infection Control Pelvic Belt ¹ – Large (frame widths 18"-22", 2-piece w/ side-release buckle)	\$95	<input type="checkbox"/>	ZCBICL	Infection Control Chest Belt ^{1,2} – Large (frame widths 18"-22", 1-piece w/ side-release buckle)	\$85		
<input type="checkbox"/>	ZPBICXL	Infection Control Pelvic Belt ¹ – XLarge (frame widths 22"-30", 2-piece w/ side-release buckle)	\$105	<input type="checkbox"/>	ZCBICXL	Infection Control Chest Belt ^{1,2} – XLarge (frame widths 22"-30", 1-piece w/ side-release buckle)	\$95		
<input type="checkbox"/>	ZPBBPS	Pelvic Belt – Bodypoint Aeromesh – Small	\$245	<input type="checkbox"/>	ZCBBPS	Chest Belt – Bodypoint Aeromesh – Small	\$255		
<input type="checkbox"/>	ZPBBP_*	Pelvic Belt – Bodypoint Aeromesh * enter size to complete code (M / L)	\$195	<input type="checkbox"/>	ZCBBP_*	Chest Belt – Bodypoint Aeromesh * enter size to complete code (M / L / XL)	\$205		
<input type="checkbox"/>	ZPBBPXL	Pelvic Belt – Bodypoint Aeromesh – XL	\$165						
Recommended sizes for Bodypoint belts	Chest belt	ZCBBPS	14"-16"W back frame	ZCBBPM	16"-18"W back frame	ZCBBPL	18"-22"W back frame	ZCBBPXL	22"-30"W back frame
	Pelvic belt	ZPBBPS	14"-16"W base frame	ZPBBPM	16"-18"W base frame	ZPBBPL	18"-22"W base frame	ZPBBPXL	22"-30"W base frame
¹ Non-absorbant belts – recommended for infection-control environments ² Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM/ZCBL/ZCBXL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.									

14. Other

<input type="checkbox"/>	Z66AS4	Anti-Tippers for 4" casters – Front (pr)	\$245	<input type="checkbox"/>	ZIHW	Whizard – winged urine deflector ³	\$67
<input type="checkbox"/>	Z66AS5	Anti-Tippers for 5" casters – Front (pr)	\$245	<input type="checkbox"/>	ZBUDBH_*	Bard Urine Drainage Bag Holder * Enter size to complete code (16" / 18" / 20" / 22")	\$155
<input type="checkbox"/>	Z66AS6	Anti-Tippers for 6" casters – Front (pr)	\$245				
<input type="checkbox"/>	Z66SPR	Anti-Tippers – telescoping for regular frame – Rear ¹ (pr)	\$110				
<input type="checkbox"/>	Z66SPLW	Anti-Tippers – telescoping for low frame – Rear (pr) for Z2LOWF4	\$110				
<input type="checkbox"/>	ZPSK1	Commode Pan Spacer Kit ² – 1" lower	\$200				
<input type="checkbox"/>	ZPSK2	Commode Pan Spacer Kit ² – 2" lower	\$215				
<input type="checkbox"/>	ZPSK3	Commode Pan Spacer Kit ² – 3" lower	\$230				
<input type="checkbox"/>	ZPSK4	Commode Pan Spacer Kit ² – 4" lower	\$245				
¹ Included with Reclining Back Option (Z515SPR/Z520SPR) ² Commode Pan Spacer Kits provide extra clearance for self-care with pan in place. Reduces toilet bowl clearance height.				³ Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place. ⁴ Must be held in place			

15. Credits and Omissions

<input type="checkbox"/> Z905	Omit Commode Pan and Mounting Brackets	-\$20	<input type="checkbox"/> Z907	Omit Foot Supports - if no other Foot / Leg Support Option is selected	-\$45
-------------------------------	--	-------	-------------------------------	--	-------

Some assembly of Raz products is required.

16. Notes:**Specifications****Raz-SP (Self Propel)**

	22" wheels	24" wheels
Weight Capacity	400 lbs	400 lbs
Overall Width	Seat width plus 7.75"	Seat width plus 7.75"
Overall Length (foot supports in middle position)	41"	42"
Seat Depth	16" or 21" depending on seat/cushion length	
Clearance Height (under seat)		
- with pan brackets	15.30"-19.30" (1" increments)	
- without pan brackets	15.93"-19.93" (1" increments)	
Clearance Width (to fit around toilet) – 18"W frame	16"	16"
Arm Support Height - above molded seat (1/2" increments)		
- 17"H Fixed Back - standard	8.5"	8.5"
- 20"H Fixed Back	9.5"	9.5"
- 15"H Adjusta-Back	5.5"-8.5"	5.5"-8.5"
- 17"H Adjusta-Back	6.5"-9.5"	6.5"-9.5"
- 20"H Adjusta-Back	7.5"-10.5"	7.5"-10.5"
- 24"H Adjusta-Back	7.5"-10.5"	7.5"-10.5"
- 20"H Reclining Adjusta-Back	7.5"-10.5"	7.5"-10.5"
- 24"H Reclining Adjusta-Back	7.5"-10.5"	7.5"-10.5"
Seat-to-Floor Height (with molded seat)	20"-24"	20"-24"
Product Weight (with 18"W seat)	47 lbs	49 lbs
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)	