



Raz-SP

MSRP – USA – USD
Effective Jan 1, 2022

Date: _____

Dealer / Location: _____ T: _____ F: _____

PO# / Client: _____ Contact: _____ Email: _____

Raz-SP – Self Propel Mobile Shower Commode Chair with height-adjustable stainless-steel frame, flip-up padded arm supports, foot supports, 22" or 24" rear wheels, front dual-locking 5" (125mm) casters and commode pan. **User weight capacity 400 lbs.** See our [Evaluation and Medical Justification Form](#).

1. Chair Codes by size – Must Choose One

Z20018	Raz-SP w/ 18"W frame – Standard base frame size, will fit over most toilets ¹	\$1,950
Z200__*	Non-Standard Width Raz-SP – Must add upcharge for custom frame width: see below * add custom width to complete code (16/20/22) – Adjusta-Back must be selected for all custom frame widths [Sec.5]	\$1,950
Z2SW16	Upcharge for custom 16"W frame ² – will not fit over most toilets	\$495
Z2SW20	Upcharge for custom 20"W frame ³	\$495
Z2SW22	Upcharge for custom 22"W frame ⁴	\$595

¹ Fits all 18"W seats
² Visco-Foam Interface seats can be modified to fit.
³ 18"W seats can be selected but must also select Seat Support Adapter Kit Z4SSKAS [Sec.2]
⁴ Z4V23N Visco Foam Interface HD seat can be modified to fit. Seat depth range (IPAS) will be 18"-20".

Custom Base Frame Length and Height

Z2LNGF3	Upcharge – custom 3" Longer Frame recommended for seat depths 20"/21", req for seat depths 21"-24"	\$395
Z2LOWF4	Upcharge for custom 4" Lower Frame – reduces seat height range to 16"-18" from standard height range of 20"-24" with molded seat – To lower chair by an additional 1" (15"-17"), add 4" caster upgrade option Z64SP [Sec.7] – Compatible only with 22" wheels [Sec.8]	\$495

2. Seats – Standard

Molded Seats (thickness ~2.5"-3")			Visco-Foam Interface Seats (thickness ~3.5")		
Z405	18"W x 16"D – Front Access	Standard	Z436BRI	18"W x 16"D – Front Bridged	\$425
Z410	18"W x 16"D – Front Bridged	N/C	Z436FRO	18"W x 16"D – Front Access	\$425
Z415	18"W x 19"D – Front Access	\$135	Z436LTO	18"W x 16"D – Left Access	\$425
Z420	18"W x 19"D – Front Bridged	\$135	Z436RTO	18"W x 16"D – Right Access	\$425
Specialty Seats			Z436RRO	18"W x 16"D – Rear Access	\$425
Z4SOM	Shower-Only Seat – Molded 18"W x 16"-21"D	\$395	Z439BRI	18"W x 19"D – Front Bridged	\$495
Z4EJRM	E&J Replica Seat – Molded – 18"W x 17"D	\$295	Z439FRO	18"W x 19"D – Front Access	\$495
Z4EJRV	E&J Replica Seat – Visco Foam – 18"W x 17"D	\$895	Z439LTO	18"W x 19"D – Left Access	\$495
Z4V23N	23"W x 18"D HD – Visco Foam Front Access – Notched at rear to maintain fore/aft adjustment range (IPAS) – 22"W frame	\$695	Z439RTO	18"W x 19"D – Right Access	\$495
			Z439RRO	18"W x 19"D – Rear Access	\$495
			Z4SSKAS	Seat Support Kit for Raz-AP/SP – for seat 2" narrower than seat frame	\$145
			Z4GSB__*	Custom Seat Bracket Spacers – for adding dump angle (pr) Specify additional height at front of seat to complete code (0.5"/1.0"/1.5")	\$295

3. Seats – Custom

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required.

Z4CS	Custom Seat (enter specs below) Select if both width and depth are ≤ 21"	\$795
Z4CEJR	Custom E&J Replica Seat (enter specs below) Select if both width and depth are ≤ 21"	\$995
Z4CSHD	Heavy Duty Custom Seats (enter specs below) Select if width and/or depth is ≥ 22"	\$1,095

Custom Specifications (for above custom seats)

Seat Frame Width:	-- 16" / 18" / 20" / 22"
Custom Seat Width:	-- 16" / 18" / 20" / 22" / 24" (cannot be more than 2" wider than frame width)
Custom Seat Depth:	-- 15" / 16" / 17" / 18" / 19" / 20" / 21" / 22" / 23" / 24" [≥ 21" requires Long Frame (Z2LNGF3)]
Access Opening:	-- FR – Front / LT – Left / RT – Right / RR – Rear / BR – Bridged
Access Opening Width:	-- 2.5" / 3.0" / 3.5" / 4.0" / 4.5" / 5.0" (3.5"W is standard)
Aperture Width:	-- 3.0"-8.0" (6.5" recommended) – enter "NA" if aperture is not required

- Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seats.
- Custom Seat width should not be narrower than frame width.
- Custom Seat can be up to 2" wider than seat frame width to provide additional thigh support without ordering wider chair frame.
- The Ischial Pelvic Alignment System (IPAS) allows 2" fore/aft adjustment of the seat (seat depth and aperture position) on all seats.
- Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, etc.

4. Fixed Back Frames (18"W only)

This is called a "Fixed" Back because the arm supports are not adjustable in height. Includes molded, vinyl back support, fixed-height flip-up padded arm supports, back angle adjustment (85°/90°/95°) with cane-style posts. Back height is measured from top of molded seat to top of back upholstery. Hip Adductors [Z7HP1/2] are not compatible with Fixed Back.

Z50818	Fixed Back – 18"W x 18"H ^{1,2}	Standard
Z5BU1818	Tension-Adjustable Back Upholstery for Fixed Back (Z508)	\$185

5. Adjusta-Back Frames


Includes height- and angle-adjustable arms with Armrest Lock Kit (Z710), tension-adjustable back upholstery, flip-up padded arm supports, back angle adjustment (85°/90°/95°). Back heights of 15"/17"/20" include cane-style posts, on which back upholstery can be positioned 2.5" higher than nominal height indicated below. e.g. 20"H can be moved up to 22.5"H. The 24"H Adjusta-Back includes straight posts, on which the upholstery cannot be raised.

Standard Width Adjust-Backs (18"W)			Non-Standard Width Adjust-Backs		
Z50918	Adjusta-Back – 15"H	\$225	* Complete Adjusta-Back order code by adding desired back frame width (16"/20"/22") Note: 18" is the standard width.		
Z51018	Adjusta-Back – 17"H	\$225	Z509__*	Adjusta-Back – 15"H non-std width	\$225
Z51518	Adjusta-Back – 20"H	\$295	Z510__*	Adjusta-Back – 17"H non-std width	\$225
Z52018	Adjusta-Back – 24"H	\$345	Z515__*	Adjusta-Back – 20"H non-std width	\$295
Z5DTL	Dual Triggers – both on left side ¹	N/C	Z520__*	Adjusta-Back – 24"H non-std width	\$345
Z5DTR	Dual Triggers – both on right side ¹	N/C	Z515SPR__*	Reclining Adjusta-Back 30° recline, 20" back height with 2-Point Reclining Arms (incl Rear Anti-tippers, but not Head Support)	\$895
Z515SPR18	Reclining Adjusta-Back 30° recline, 20" back height with 2-Point Reclining Arms (incl Rear Anti-tippers, but not Head Support)	\$895	Z520SPR__*	Reclining Adjusta-Back 30° recline, 24" back height with 2-Point Reclining Arms (incl Rear Anti-tippers but not Head Support)	\$995
Z520SPR18	Reclining Adjusta-Back 30° recline, 24" back height with 2-Point Reclining Arms (incl Rear Anti-tippers but not Head Support)	\$995			

¹ Dual Tilt Triggers option positions both triggers on one side, allowing for single-handed activation. Available for all back widths.

6. Custom Back Frame Width Narrower Than Chair Frame Width – Must choose back frame [Sec.5] (configurations below are not available for recliner backs)

NOTE: For back frame width variances on recliner frames, please request custom quotation.

Z5NB__ __ *	Custom Back Frame 2" Narrower than base frame	\$345	Base Frame Width	16" ¹	18" ²	20" ³	22" ⁴	
			Back Frame Width	14"	16"	18"	20"	
<p>* Enter base and back frame widths to complete code. eg. Z5NB2018 specifies a 20"W chair frame with an 18"W back frame</p>			<p>Consider also using Arm Spacer Kit (Z718/Z719) as another solution for providing more width between the arms.</p>					
<p>¹ Requires Custom Seat</p> <p>² Compatible with all 18"W seats but with reduced seat depth range. Also compatible with Custom or modified Visco-Foam Interface seats.</p> <p>³ Custom Seat recommended. Also compatible with all 18"W seats with Seat Support Kit Z4SSKAS [Sec.2].</p> <p>⁴ Custom Seat recommended. Also compatible with modified Z4V23N (18" seat depth with limited adjustability)</p>								

7. Back Support Options – These back supports replace the standard, tension-adjustable upholstery and are recommended for infection-control environments.

Z5S14	14"W Symphony Back ¹ – recessed solid back	\$395	Z749 Harmony Back ² (contoured solid back – 18"W x 13"H) Fits 17" or 20" back heights, 18"W back widths only	\$245
Z5S16	16"W Symphony Back ¹ – recessed solid back	\$395		
Z5S18	18"W Symphony Back ¹ – recessed solid back	\$395		
Z5S20	20"W Symphony Back ¹ – recessed solid back	\$495		
Z5S22	22"W Symphony Back ¹ – recessed solid back	\$595		
			<p>¹ Not compatible w/ 15" back height – Symphony Back panels are 13"H</p> <p>² Not compatible w/ 15" and 24" back heights. Not recommended for reclining backs.</p>	

8. Rear Wheel, Caster Options and Anti-Tippers

Z612R	22" Diameter Wheels w/ Textured Handrims	Std	Z660	Quick-Release Axles (pr)	\$95
Z614R	24" Diameter Wheels w/ Textured Handrims	N/C	Z674	User-Actuated Front Directional Caster Locks w/ extension handles (pr)	\$345
Z6WLE	Wheel Lock with Extensions (pr)	\$95	Z64SP	4" Caster Upgrade affects STF range (2x)	\$140
Z66SPR	Rear Anti-Tippers for regular frame ¹ (pr)	\$110	Z622VHR	Vinyl-Coated Handrims for 22" (pr)	\$265
Z66SPLW	Rear Anti-Tippers for low frame (pr) Z2LOWF4	\$110	Z624VHR	Vinyl-Coated Handrims for 24" (pr)	\$265
Z66AS4	Front Anti-Tippers for 4" casters (pr)	\$245	Z622PHR	Vertical Projection Handrims – 22" (pr)	\$465
Z66AS5	Front Anti-Tippers for 5" casters (pr)	\$245	Z624PHR	Vertical Projection Handrims – 24" (pr)	\$495
Z66AS6	Front Anti-Tippers for 6" casters (pr)	\$245			

9. Arm Supports and Accessories

Z5BAR13	Flip-Up Arm Supports – 13”L pad (pr)	Std	Z7PAML	Pivoting Arm Mount ³ – Left (for Z7FL/Z7TR)	\$165
Z5BAR10	Flip-Up Arm Supports – 10”L pad (pr) (not compatible with Z7FS/Z7FL/Z7TR / Z7PAM)	N/C	Z7PAMR	Pivoting Arm Mount ³ – Right (for Z7FL/Z7TR)	\$165
Z710	Arm Support Lock Kit (for Z505/Z506) (pr)	\$75	Z718	Arm Support Spacer Kit – 1/2” per side (for Z509/Z510/Z515/Z520) (pr)	\$60
Z52PAS	2-Point Flip-up Arm Supports ¹ (pr) (standard with reclining backs)	\$345	Z719	Arm Support Spacer Kit – 1” per side (for Z509/Z510/Z515/Z520) (pr)	\$115
Z7FLL	Long Molded Flat Armpad ² – Left (14”L x 3.25”W, fits Z5BAR13)	\$95	Z7PHG	Pivoting Hand Grips (pr)	\$395
Z7FLR	Long Molded Flat Armpad ² – Right (14”L x 3.25”W, fits Z5BAR13)	\$95	Z7APS_ _*	Anterior Posture Support Bar * Enter back frame width to complete code (14”/16”/18”) ^{3, 4, 5}	\$195
Z7FSL	Short Molded Flat Armpad – Left (10”L x 3”W, fits Z5BAR10)	\$95	Z7APS_ _*	Anterior Posture Support Bar * Enter back frame width to complete code (20”/22”) ^{3, 4, 5}	\$225
Z7FSR	Short Molded Flat Armpad – Right (10”L x 3”W, fits Z5BAR10)	\$95			
Z7TRL	Molded Arm Trough ² – Left (18.5”L x 4.5”W, fits Z5BAR13)	\$110			¹ For chairs with same back width and frame width only
Z7TRR	Molded Arm Trough ² – Right (18.5”L x 4.5”W, fits Z5BAR13)	\$110			² May require Z7PAM for arm support to flip up and/or for lateral supports to swing away
					³ Not compatible with Z5BAR10
					⁴ Not compatible with Z7FL/Z7FS/Z7TR
					⁵ Special length required with Z718 (+1”) and Z719 (+2”)

10. Foot and Leg Supports – Standard Foot Supports (Z705N) are included with the chair. The selection of any Foot / Leg Support Option below automatically replaces the relevant Standard Foot Support. e.g. By selecting Elevating Leg Support – Left (Z7ELRL), the chair will be supplied with one Foot Support (Z705N) on the right side and a Z7ELRL on the left side.

Z705N	Foot Supports 16”-20.5” (pr)	Std	Z7ELRL	Elevating Leg Support ³ – Left (includes Z7ADM)	\$475
Z705N16	Foot Supports – 16”W frames 16”-20.5” (pr)	N/C	Z7ELRR	Elevating Leg Support ³ – Right (includes Z7ADM)	\$475
Z7FET2	Footrest Extension Tubes 2” longer ¹ (pr)	\$50	Z7ELROFL	Forward Offset Elevating Leg Support ⁴ – Left (includes Z7ADM)	\$525
Z7MFX8	MFX8 Foot Supports 8”-12” (pr)	\$75	Z7ELROFR	Forward Offset Elevating Leg Support ⁴ – Right (includes Z7ADM)	\$525
Z7MFX12	MFX12 Foot Supports 12”-16.5” (pr)	\$175	Z7RLSL	Universal Residual Limb Support – Left	\$395
Z7ADM	Adj Angle/Depth Molded Footplates ² (pr)	\$190	Z7RLSR	Universal Residual Limb Support – Right	\$395
Z7FET2A	Footrest Extension Tubes for Z7ADM (pr) (2” longer)	\$50			
Z7ADABS	Angle/Depth Adj ABS Footplates (pr) for 16”W frames	\$190			
Z7VS18	V-Style Foot Support (1-piece) – for 18”W base frames	\$195			
Z7VS_ _*	V-Style Foot Support (1-piece) * Enter frame width to complete code (16”/18”/20”/22”)	\$295			
Z7FBAS	Flip-Back Foot Support (1-piece)	\$195			
					³ Z7ELRL/R compatible with frame widths ≥ 18” and adjusted seat depths ≤ 18” or ≤ 21” w/ 3” longer seat frames (Z2LNGF3)
					⁴ Z7ELROFL/R compatible with frame widths ≥ 18” for adjusted seat depths 19”-20” on standard frame or 22”-23” w/ 3” longer seat frames (Z2LNGF3)
					¹ Not compatible with Z7MFX8/Z7ADM/Z7VS/Z7FBAS
					² Not compatible with Z20016 or Z7MFX8

11. Foot and Leg Support Accessories

Z7FSLK	Foot Support Lock (pr)	\$95	Z7NCS	Neoprene Covers for Standard Footplates (pr)	\$90
Z7CS__*	Snap-On Calf Strap – Bodypoint Aeromesh* Enter frame width to complete code (16"/18"/20"/22")	\$95	Z7NC16	Neoprene Covers for narrow molded footplates (pr)	\$90
Z7CSIC__*	Snap-On Infection Control Calf Strap * Enter seat frame width to complete code (16"/18"/20")	\$125	Z7NCA	Neoprene Covers for Adj Angle/Depth Molded Footplates (pr)	\$90
Z7CSIC22	Snap-On Infection Control Calf Strap – 22"W frames	\$140	Z7NCABS	Neoprene Covers for ABS Adj Angle/Depth Footplates Z7ADABS (pr)	\$90
Z737CPL	Calf Panel – Bodypoint Aeromesh ¹ Fits frame widths (16"/18"/20")	\$175	Z7OFLM	Lateral Offset Foot Support Receiver – moves foot support 2" laterally or medially Quantity (1 or 2) __	\$160
Z706	H-Strap – for 18"/20"W frames – hook and loop attachment	\$90			
Z712	H-Strap – for 22"W frames – hook and loop attachment	\$110			
Z708ADJ	Adjustable Heel Loops (pr) – not available for Z7VS / Z7FBAS	\$60			

¹ Not compatible with Z7MFX8 / Z7FBAS / Z7VS

12. Positioning Pads

Z7HSASS	Head Support w/ Small Pad (10" x 4")	\$495	Z7LGAS1L	Leg Adductor w/ Small Pad ¹ – Left	\$320
Z7HSASL	Head Support w/ Large Pad (12" x 6")	\$545	Z7LGAS1R	Leg Adductor w/ Small Pad ¹ – Right	\$320
Z7OHIP	Lateral Offset Head Support Interface Plate– up to 22° of adjustment	\$95	Z7LGAS2L	Leg Adductor w/ Large Pad ¹ – Left	\$320
Z7LTS1L	Lateral Thoracic Support w/ Small Pad ^{1,2} – Left	\$320	Z7LGAS2R	Leg Adductor w/ Large Pad ¹ – Right	\$320
Z7LTS1R	Lateral Thoracic Support w/ Small Pad ^{1,2} – Right	\$320	Z7HP1L	Hip Adductor w/ Small Pad ^{1,3} – Left	\$320
Z7LTS2L	Lateral Thoracic Support w/ Large Pad ^{1,2} – Left	\$320	Z7HP1R	Hip Adductor w/ Small Pad ^{1,3} – Right	\$320
Z7LTS2R	Lateral Thoracic Support w/ Large Pad ^{1,2} – Right	\$320	Z7HP2L	Hip Adductor w/ Large Pad ^{1,3} – Left	\$320
Z742L	Lateral/Adductor Extension Plate (5.5" x 3.5") – Left	\$25	Z7HP2R	Hip Adductor w/ Large Pad ^{1,3} – Right	\$320
Z742R	Lateral/Adductor Extension Plate (5.5" x 3.5") – Right	\$25			

¹ Adjustable Pad Sizes: Small – 5.5"x 3.5" / Large – 6.75" x 4.3".

² Pads on Lateral Thoracic Supports can be adjusted medially up to 3" per side. For additional medial offset, specify Z742 extension plate.

³ Fits Adjusta-Back only

13. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts				Chest Belts					
ZPBS	Pelvic Belt – Small (frame widths 14”-16”, 2-piece w/ side-release buckle)	\$35		ZCBS	Chest Belt ² – Small (frame widths 14”-16”, 1-piece w/ hook and loop closure)	\$35			
ZPBM	Pelvic Belt – Medium (frame widths 16”-18”, 2-piece w/ side-release buckle)	\$35		ZCBM	Chest Belt ² – Medium (frame widths 16”-18”, 1-piece w/ hook and loop closure)	\$35			
ZPBL	Pelvic Belt – Large (frame widths 18”-22”, 2-piece w/ side-release buckle)	\$45		ZCBL	Chest Belt ² – Large (frame widths 18”-22”, 1-piece w/ hook and loop closure)	\$45			
ZPBXL	Pelvic Belt – XLarge (frame widths 22”-30”, 2-piece w/ side-release buckle)	\$55		ZCBXL	Chest Belt ² – XLarge (frame widths 22”-30”, 1-piece w/ hook and loop closure)	\$55			
ZPBICS	Infection Control Pelvic Belt ¹ – Small (frame widths 14”-16”, 2-piece w/ side-release buckle)	\$75		ZCBICS	Infection Control Chest Belt ^{1,2} – Small (frame widths 14”-16”, 1-piece w/ side-release buckle)	\$75			
ZPBICM	Infection Control Pelvic Belt ¹ – Medium (frame widths 16”-18”, 2-piece w/ side-release buckle)	\$85		ZCBICM	Infection Control Chest Belt ^{1,2} – Medium (frame widths 16”-18”, 1-piece w/ side-release buckle)	\$75			
ZPBICL	Infection Control Pelvic Belt ¹ – Large (frame widths 18”-22”, 2-piece w/ side-release buckle)	\$95		ZCBICL	Infection Control Chest Belt ^{1,2} – Large (frame widths 18”-22”, 1-piece w/ side-release buckle)	\$85			
ZPBICXL	Infection Control Pelvic Belt ¹ – XLarge (frame widths 22”-30”, 2-piece w/ side-release buckle)	\$105		ZCBICXL	Infection Control Chest Belt ^{1,2} – XLarge (frame widths 22”-30”, 1-piece w/ side-release buckle)	\$95			
ZPBBPS	Pelvic Belt – Bodypoint Aeromesh – Small	\$245		ZCBBPS	Chest Belt – Bodypoint Aeromesh – Small	\$255			
ZPBBP_*	Pelvic Belt – Bodypoint Aeromesh * enter size to complete code (M / L)	\$195		ZCBBP_*	Chest Belt – Bodypoint Aeromesh * enter size to complete code (M / L / XL)	\$205			
ZPBBPXL	Pelvic Belt – Bodypoint Aeromesh – XL	\$165							
Recommended sizes for Bodypoint belts	Pelvic belt	ZPBBPS	14”/16”W base frame	ZPBBPM	16”/18”W base frame	ZPBBPL	18”-22”W base frame	ZPBBPXL	22”-26”W base frame
	Chest belt	ZCBBPS	14”/16”W back frame	ZCBBPM	16”/18”W back frame	ZCBBPL	18”/20”W back frame	ZCBBPXL	22”-26”W back frame
¹ Non-absorbent belts – recommended for infection-control environments ² Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM/ZCBL/ZCBXL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.									

14. Other

ZPSK1	Commode Pan Spacer Kit ² – 1” lower	\$200	ZIHW	Whizard – winged urine deflector ³	\$67
ZPSK2	Commode Pan Spacer Kit ² – 2” lower	\$215	ZBUDBH_*	Bard Urine Drainage Bag Holder * Enter size to complete code (16” / 18” / 20” / 22”)	\$155
ZPSK3	Commode Pan Spacer Kit ² – 3” lower	\$230			
ZPSK4	Commode Pan Spacer Kit ² – 4” lower	\$245			
¹ Included with Reclining Back Option (Z515SPR/Z520SPR) ² Commode Pan Spacer Kits provide extra clearance for self-care with pan in place. Reduces toilet bowl clearance height.			³ Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place. ⁴ Must be held in place		

15. Credits and Omissions

Z905	Omit Commode Pan and Mounting Brackets	-\$20	Z907	Omit Foot Supports - if no other Foot / Leg Support Option is selected	-\$45
------	--	-------	------	--	-------

Some assembly of Raz products is required.

16. Notes:**Specifications****Raz-SP (Self Propel)**

	22" wheels	24" wheels
Weight Capacity	400 lbs	400 lbs
Overall Width	Seat width plus 7.75"	Seat width plus 7.75"
Overall Length (foot supports in middle position)	41"	42"
Seat Depth	16" or 21" depending on seat/cushion length	
Clearance Height (under seat)		
- with pan brackets	15.25"-19.25" (1" increments)	
- without pan brackets	16"-20" (1" increments)	
Clearance Width (to fit around toilet) – 18"W frame	16"	16"
Arm Support Height - above molded seat (1/2" increments)		
- 18"H Fixed Back - standard	8.5"	8.5"
- 15"H Adjusta-Back	5.5"-8.5"	5.5"-8.5"
- 17"H Adjusta-Back	6.5"-9.5"	6.5"-9.5"
- 20"H Adjusta-Back	7.5"-10.5"	7.5"-10.5"
- 24"H Adjusta-Back	7.5"-10.5"	7.5"-10.5"
- 20"H Reclining Adjusta-Back	7.5"-10.5"	7.5"-10.5"
- 24"H Reclining Adjusta-Back	7.5"-10.5"	7.5"-10.5"
Seat-to-Floor Height (with molded seat)	20"-24"	20"-24"
Product Weight (with 18"W seat)	47 lbs	49 lbs
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)	