

Accessibility for everyone...

TM

The Convert-able Table



The Convert-able Table is so adaptable that it can become a desk, an art easel or even a sensory stim bowl. It is appropriate for adults as well as children. The push button height adjustment allows the Convert-Able table to be used with wheelchairs and other adaptive equipment including standing frames, seating system and gait trainers.

The table has three interchangeable tops that simply snap on and off. The Desk Top, made of durable melamine, is suitable in the classroom, resource center or home setting. The large work surface, with up to 90 degrees of angled tilt, makes the Convert-Able Tale appropriate for therapy exercises as well.

Whether using watercolors, finger paints, crayons or clay, both skilled and budding artists will find no limit to the kinds of creative activities the Art Easel Top can accommodate.

The durable scratch sistant surface wipes clean. The ample work surface keeps paper unwrinkled on a roll while horizontal guides hold the paper in place, even under vigorous brush work. The accessory art bin, which attaches to the table, will help keep art supplies organized.

Encourage imaginative play, watch goldfish swim, make mud pies or play with soap bubbles in the Sensory Stim Bowl. Made of clear acrylic, the bowl is large enough for play or therapy activities,

No other Table but the Convert-able by R.E.A.L. Design offers so much versatility!



The Convert-able Table

TM



Can be used by all - including elementary students, seniors and others using wheelchairs or other adaptive seating

Height adjustable with the push of a button

The work surface adjusts in angle from 0 to 90 degrees using a single lever

Ample size desk top can be used as a computer work station

Compact, portable base is easily movable from place to place

Once in place, dual lock casters position table securely

Custom school colors are available on the table frame

Accessory items available include book bin, pencil holder and the Big Bend Clip (as pictured on other side)

	Inside width of frame	26"	66cm.
	Convert-Able Table Height	21" - 37"	53 - 94 cm.
	Desk Top Dimension	27" X 31"	68.5 X 79 cm.
	Art Easel Dimensions	27" x 31"	68.5 X 79 cm.
	Table Weight w/o Tops	45 lbs.	20.5 kg
	Wheel Diameter	3"	7.6 cm.
	Outside Frame Dimensions	29"x 30"	74x76 cm

