

# REBOTEC

## AUGSBURG

### Shower commode chair for children



- Custom design for children for use in the shower, over the toilet and at bedside.
- Soft seat and backrest
- Upholstered seat with hygiene opening
- Height adjustable, swinging and detachable foot support
- Swinging armrests for easy side transfer

#### MPN

Augsburg 4" Castors  
339.04.97

Augsburg height adjustable head-rest  
444.03.39

Augsburg seat width and depth adjustment kit  
444.00.97

Augsburg 24" Self Propelled  
339.24.97

#### ACCESSORIES

Bucket Pail  
301.00

Safety Belt  
322.30.10

IV Pole  
303.62

Armrest Cushion  
440.10.32

CHAIR DEPTH, WIDTH 87x53 cm

SEAT DEPTH, WIDTH, HEIGHT 35x40x54 cm

TOILET HEIGHT WITHOUT BUCKET DEVICE 47 cm

TOILET HEIGHT WITH BUCKET DEVICE 43 cm

TOILET WIDTH 33 cm

MAX BODY WEIGHT 80 kg



Australian Government  
Department of Health  
Therapeutic Goods Administration

ARTG 297301

