

2021 order form for large

Use dimension chart to select appropriate size.



		Standard R860 large	Hi/lo R870 large
* Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R861 \$660 <input type="checkbox"/> R871 \$1,260 <input type="checkbox"/>
		Adjustable w/ back spring	R862 \$765 <input type="checkbox"/> R872 \$1,365 <input type="checkbox"/>
* Base	CHOOSE	Without spring	R873 \$405 <input type="checkbox"/> R883 \$2,730 <input type="checkbox"/>
		With spring	R874 \$485 <input type="checkbox"/> N/A
* Arm supports (pair)	CHOOSE	Armrests	R868 \$235 <input type="checkbox"/> R868 \$235 <input type="checkbox"/>
		Forearm prompts	R879 \$610 <input type="checkbox"/> R879 \$610 <input type="checkbox"/>
* Adjustable legs	CHOOSE	Short	R808 \$65 <input type="checkbox"/> N/A
		Long	R894 \$75 <input type="checkbox"/> N/A
		Short with casters	R895 \$285 <input type="checkbox"/> N/A
		Long with casters	R896 \$310 <input type="checkbox"/> N/A
* Pads	CHOOSE	Tan	R875 \$245 <input type="checkbox"/> R875 \$245 <input type="checkbox"/>
		Blue	R876 \$245 <input type="checkbox"/> R876 \$245 <input type="checkbox"/>
		Red	R877 \$245 <input type="checkbox"/> R877 \$245 <input type="checkbox"/>
		Pink	R977 \$245 <input type="checkbox"/> R977 \$245 <input type="checkbox"/>
		Purple	R978 \$245 <input type="checkbox"/> R978 \$245 <input type="checkbox"/>
		Green	R979 \$245 <input type="checkbox"/> R979 \$245 <input type="checkbox"/>
Headrest	CHOOSE	Contoured	R897 \$230 <input type="checkbox"/> R897 \$230 <input type="checkbox"/>
		Flat	R898 \$215 <input type="checkbox"/> R898 \$215 <input type="checkbox"/>
		Adjustable winged	R881 \$285 <input type="checkbox"/> R881 \$285 <input type="checkbox"/>
Laterals (pair)	CHOOSE	Small	R828 \$190 <input type="checkbox"/> R828 \$190 <input type="checkbox"/>
		Large	R838 \$205 <input type="checkbox"/> R838 \$205 <input type="checkbox"/>
		Large with strap	R839 \$245 <input type="checkbox"/> R839 \$245 <input type="checkbox"/>
Push handles (recommended for all chairs with wheels)		R878 \$145 <input type="checkbox"/>	Included
Tray (requires armrests)		R669 \$275 <input type="checkbox"/>	R669 \$275 <input type="checkbox"/>
Handhold		R649 \$25 <input type="checkbox"/>	R649 \$25 <input type="checkbox"/>
Chest strap (wide)		R807 \$90 <input type="checkbox"/>	R807 \$90 <input type="checkbox"/>
Butterfly harness	CHOOSE	Regular	R889 \$140 <input type="checkbox"/> R889 \$140 <input type="checkbox"/>
		Slim cut	R826 \$140 <input type="checkbox"/> R826 \$140 <input type="checkbox"/>
Thigh belt		R813 \$90 <input type="checkbox"/>	R813 \$90 <input type="checkbox"/>
Pelvic harness		R846 \$170 <input type="checkbox"/>	R846 \$170 <input type="checkbox"/>
Hip guides (pair)		R867 \$200 <input type="checkbox"/>	R867 \$200 <input type="checkbox"/>
Abductor	CHOOSE	Small	R884 \$120 <input type="checkbox"/> R884 \$120 <input type="checkbox"/>
		Large	R891 \$140 <input type="checkbox"/> R891 \$140 <input type="checkbox"/>
Adductors (pair)		R892 \$190 <input type="checkbox"/>	R892 \$190 <input type="checkbox"/>
Leg prompt		R859 \$225 <input type="checkbox"/>	R859 \$225 <input type="checkbox"/>
Footboard (recommended for chairs with long legs)		R857 \$470 <input type="checkbox"/>	Included
Ankle straps (pair)		R855 \$90 <input type="checkbox"/>	R855 \$90 <input type="checkbox"/>
Sandals (pair)		R964 \$105 <input type="checkbox"/>	R964 \$105 <input type="checkbox"/>
Wedges (pair, require sandals)		R699 \$70 <input type="checkbox"/>	R699 \$70 <input type="checkbox"/>
Lumbar and seat support kit		R888 \$85 <input type="checkbox"/>	R888 \$85 <input type="checkbox"/>
Backrest filler pad		R856 \$120 <input type="checkbox"/>	R856 \$120 <input type="checkbox"/>
Whitmyer headrest adapter		R804 \$55 <input type="checkbox"/>	R804 \$55 <input type="checkbox"/>
Total		R860 \$ _____	R870 \$ _____

* Required

FREE SHIPPING •

Room to grow

Large chair at its smallest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 9½"
- 3 Seat depth 15"
- 4 Seat width 11"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 11¼"




Large chair at its largest dimensions

- 1 Backrest height 24"
- 2 Distance between laterals 14"
- 3 Seat depth 20"
- 4 Seat width 17"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 21¼"




Large standard base
R860 Rifton Activity Chair

User dimensions (inches)	
Height	50–74
 <p>Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p>	
Item dimensions (inches)	
	R860 large Standard base
Frame width	short legs: 26¾ long legs: 26¾
	short legs w/ casters: 27¾ long legs w/ casters: 27¾
Seat height above floor	short legs: 16½–20½ long legs: 20–24
	short legs w/ casters: 16½–20½ long legs w/ casters: 20–24
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	11¼–21¼
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lb)	250



Large hi/lo base
R870 Rifton Activity Chair

User dimensions (inches)	
Height	50–74
 <p>Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p>	
Item dimensions (inches)	
	R870 large Hi/lo base
Frame width	29¾
Seat height above floor	16–24
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	11¼–21¼
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lb)	250