

2024 order form


Use dimension chart to select appropriate size.


			S410 size 1			S420 size 2				
* Frame			S411	\$1,900	<input type="checkbox"/>	S421	\$2,500	<input type="checkbox"/>		
* Color			CHOOSE	Blue	<input type="checkbox"/>	Blue	<input type="checkbox"/>			
				Red	<input type="checkbox"/>	Red	<input type="checkbox"/>			
* CHOOSE CONFIGURATION	Prone (with laterals, pelvic harness) HCPCS code: E0638	Prone knee cuffs (pair)	S412	\$500	<input type="checkbox"/>	S422	\$510	<input type="checkbox"/>		
	Supine (with upper trunk support, laterals) HCPCS code: E0638	Head laterals (pair)	S413	\$550	<input type="checkbox"/>	S423	\$560	<input type="checkbox"/>		
		Butterfly harness	S497	\$200	<input type="checkbox"/>	S497	\$200	<input type="checkbox"/>		
		* Hip strap	CHOOSE	Standard	S456	\$300	<input type="checkbox"/>	S457	\$310	<input type="checkbox"/>
				Rotation control	S472	\$300	<input type="checkbox"/>	S473	\$310	<input type="checkbox"/>
				Rotation control	S476	\$300	<input type="checkbox"/>	S477	\$310	<input type="checkbox"/>
			* Soft knee straps (pair)	S401	\$200	<input type="checkbox"/>	S402	\$210	<input type="checkbox"/>	
	Multi-Position (with upper trunk support, laterals, pelvic harness) HCPCS code: E0641			S414	\$850	<input type="checkbox"/>	S424	\$870	<input type="checkbox"/>	
		Head laterals	S497	\$200	<input type="checkbox"/>	S497	\$200	<input type="checkbox"/>		
		Butterfly harness	S456	\$300	<input type="checkbox"/>	S457	\$310	<input type="checkbox"/>		
		* Hip strap	CHOOSE	Standard	S472	\$300	<input type="checkbox"/>	S473	\$310	<input type="checkbox"/>
				Rotation control	S476	\$300	<input type="checkbox"/>	S477	\$310	<input type="checkbox"/>
			* Soft knee straps (pair)	S401	\$200	<input type="checkbox"/>	S402	\$210	<input type="checkbox"/>	
			Prone knee cuffs (pair)	S485	\$200	<input type="checkbox"/>	S485	\$200	<input type="checkbox"/>	
* Trunk strap	CHOOSE	Standard	S462	\$300	<input type="checkbox"/>	S463	\$310	<input type="checkbox"/>		
		Rotation control	S466	\$300	<input type="checkbox"/>	S467	\$310	<input type="checkbox"/>		
* Leg type	CHOOSE	Standard	S415	\$650	<input type="checkbox"/>	S425	\$670	<input type="checkbox"/>		
		Advanced	N/A			S426	\$970	<input type="checkbox"/>		
Tray			S491	\$550	<input type="checkbox"/>	S491	\$550	<input type="checkbox"/>		
Additional trunk laterals (pair)			S494	\$200	<input type="checkbox"/>	S494	\$200	<input type="checkbox"/>		
Medial/lateral thigh supports (pair)			S482	\$200	<input type="checkbox"/>	S482	\$200	<input type="checkbox"/>		
Additional medial/lateral thigh supports (pair)			S482	\$200	<input type="checkbox"/>	S482	\$200	<input type="checkbox"/>		

*** Required**

The HCPCS codes indicated above are suggestions only, based on knowledge of our products and the HCPCS definitions established by CMS. It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. Actual coding and coverage varies by payer; some may accept or require alternative HCPCS codes, including miscellaneous codes, to ensure access for their beneficiaries.

Stander dimensions

User dimensions (inches)	S410 Size 1	S420 Size 2
Height	28–42	38–53
 <p>Key user dimension: height Select the appropriate stander by the user's overall height. Choose the model that allows for growth.</p> <p>Important: User's weight must not exceed the maximum working load.</p>		
Item dimensions (inches)	S410 Size 1	S420 Size 2
Base length x width	32½ x 21½	39 x 27
Height of horizontal body support above floor	22–28	22½–36½
Sandal to knee height	6¾–11	10–15¼
Knee to hip height	6¼–10½	8½–13½
Sandal to top of headrest	31½–43	41–55
Width between trunk laterals	5½–10	7–11½
Width between medial/lateral thigh supports	2½–4	3¼–5
Basic item weight PRONE (lb)	47	72
Basic item weight SUPINE (lb)	52	78
Max. working load (lb)	75	100



Size 1 range of adjustment

Size 2 range of adjustment