


How to buy

- 1. Measure** your client to determine the size you need.
(See dimension chart for help.)
- 2. Select** the appropriate size.
- 3. Choose** the accessories you need.
- 4. Add** up the total cost.

User dimensions (inches)	E420 small	E430 large
Height	30–50	46–72
 Key user dimension: height Select the appropriate stander by the user's overall height. Choose the model that allows for growth.		
	Important: User's weight must not exceed the maximum working load.	
Item dimensions (inches)	E420 small	E430 large
Board length and width	53 x 16½	73 x 19
Distance between laterals	7–13	9–16
Height when horizontal	28½	21–30
Height when vertical	54½	74
Base length and width	37½ x 27	46½ x 29½
Max. working load (lbs)	100	250



Supine Stander includes padding, footboard, adjustable armrests, hip strap, chest strap, leg straps, one pair of laterals, one pair of head blocks and double-locking casters.

		E420 small	E430 large	
		E40	\$2,600 <input type="checkbox"/>	
		E41	\$4,100 <input type="checkbox"/>	
Optional accessories	Additional laterals (pair)	Small	E504 \$135 <input type="checkbox"/> N/A	
		Large	N/A E439 \$195 <input type="checkbox"/>	
	Tray	E404	\$260 <input type="checkbox"/> E414 \$320 <input type="checkbox"/>	
	Adjustable abduction wedge	E945	\$215 <input type="checkbox"/> E945 \$215 <input type="checkbox"/>	
		Round abductor	E933	\$77 <input type="checkbox"/> E933 \$77 <input type="checkbox"/>
	Collar	E934	\$64 <input type="checkbox"/> E934 \$64 <input type="checkbox"/>	
	Sandals (pair)	Small	E422 \$82 <input type="checkbox"/> N/A	
		Medium	E433 \$90 <input type="checkbox"/> E433 \$90 <input type="checkbox"/>	
		Large	N/A E432 \$93 <input type="checkbox"/>	
	Wedges (pair, requires sandals)	Small	E424 \$59 <input type="checkbox"/> N/A	
Medium		E435 \$64 <input type="checkbox"/> E435 \$64 <input type="checkbox"/>		
Large		N/A E434 \$67 <input type="checkbox"/>		
Hand Anchor (For more hand/arm support options see pp. 98–101.)	K821	\$140 <input type="checkbox"/> K821 \$140 <input type="checkbox"/>		
Total	E420	\$ _____	E430	\$ _____

FREE SHIPPING