


2021 order form

Use dimension chart to select appropriate size.



User dimensions (inches)	R120 small	R130 medium	R140 large
Inside leg length	17 – 22	22 – 28	25 – 35
 Key user dimension: inside leg length The extended leg should reach from seat to pedal comfortably when both feet and torso are strapped into place. Important: User's weight must not exceed the maximum working load. Rifton tricycles are not intended for clients with strong reflexes or poorly developed trunk balance and head control. Consult your client's therapist before ordering.			
	Item dimensions (inches)	R120 small	R130 medium
Overall height (at lowest position)	28½	31	36½
Overall length	46	52	65
Overall width	27½	27½	31½
Wheel diameter	13	15	19
Seat center to extended pedal	17 – 22	22 – 28	25 – 35
Seat to top of trunk support	15½ – 22	15½ – 22	15½ – 22
Gear ratio (approximate)	1:1	1:1	1:1
Item weight (lb)	47	56	65
Distance between trunk laterals	8 – 14	8 – 14	8 – 14
Max. working load (lb)	150	160	200

		R120 small	R130 medium	R140 large
* Tricycle frame includes self-leveling pedals with straps, backrest tube with handle, tote box for storage and handbrake for the large only. Specify color: Red, blue, pink or lime.	CHOOSE	R112 \$1,030 <input type="checkbox"/>	R113 \$1,170 <input type="checkbox"/>	R114 \$1,980 <input type="checkbox"/>
		Red <input type="checkbox"/>	Red <input type="checkbox"/>	Red <input type="checkbox"/>
		Blue <input type="checkbox"/>	Blue <input type="checkbox"/>	Blue <input type="checkbox"/>
		Pink <input type="checkbox"/>	Pink <input type="checkbox"/>	Pink <input type="checkbox"/>
		Lime <input type="checkbox"/>	Lime <input type="checkbox"/>	Lime <input type="checkbox"/>
* Seat (includes seatbelt)	CHOOSE	Small R122 \$210 <input type="checkbox"/>	R132 \$210 <input type="checkbox"/>	N/A
	Large R128 \$355 <input type="checkbox"/>	R138 \$360 <input type="checkbox"/>	(included)	
* Handlebar	CHOOSE	Straight R126 \$135 <input type="checkbox"/>	R136 \$145 <input type="checkbox"/>	R146 \$225 <input type="checkbox"/>
	Loop R129 \$135 <input type="checkbox"/>	R139 \$145 <input type="checkbox"/>	R149 \$225 <input type="checkbox"/>	
* Back support	CHOOSE	Trunk support system R147 \$485 <input type="checkbox"/>	R147 \$485 <input type="checkbox"/>	R147 \$485 <input type="checkbox"/>
	Backrest pad R145 \$110 <input type="checkbox"/>	R145 \$110 <input type="checkbox"/>	R145 \$110 <input type="checkbox"/>	
Headrest (requires trunk support system)	CHOOSE	Winged R159 \$230 <input type="checkbox"/>	R159 \$230 <input type="checkbox"/>	R159 \$230 <input type="checkbox"/>
	Contoured R156 \$220 <input type="checkbox"/>	R156 \$220 <input type="checkbox"/>	R156 \$220 <input type="checkbox"/>	
	Flat R157 \$210 <input type="checkbox"/>	R157 \$210 <input type="checkbox"/>	R157 \$210 <input type="checkbox"/>	
Handbrakes (factory installation recommended)		N/A	R125 \$110 <input type="checkbox"/>	(included)
Communication tray	CHOOSE	Small K225 \$190 <input type="checkbox"/>	K225 \$190 <input type="checkbox"/>	K225 \$190 <input type="checkbox"/>
	Large N/A	K245 \$195 <input type="checkbox"/>	K245 \$195 <input type="checkbox"/>	
Abductor		R121 \$90 <input type="checkbox"/>	R131 \$150 <input type="checkbox"/>	R141 \$170 <input type="checkbox"/>
Front guide bar		R135 \$135 <input type="checkbox"/>	R135 \$135 <input type="checkbox"/>	R135 \$135 <input type="checkbox"/>
Rear steering bar		R123 \$235 <input type="checkbox"/>	R123 \$235 <input type="checkbox"/>	R143 \$260 <input type="checkbox"/>
Front pulley		R144 \$45 <input type="checkbox"/>	R144 \$45 <input type="checkbox"/>	R144 \$45 <input type="checkbox"/>
Additional backrest pad		R145 \$110 <input type="checkbox"/>	R145 \$110 <input type="checkbox"/>	R145 \$110 <input type="checkbox"/>
Stationary stand		R155 \$120 <input type="checkbox"/>	R155 \$120 <input type="checkbox"/>	R155 \$120 <input type="checkbox"/>
Total		R120 \$ _____	R130 \$ _____	R140 \$ _____

* Required

FREE SHIPPING