

# The Support Station from Rifton

## See what a difference 90° can make

### Upright. Dignified. Independent. Finally!

The Support Station from Rifton is a wall-mounted innovation that enables individuals with disabilities to support themselves upright during diapering or personal care. It brings new dignity, efficiency and therapeutic benefits to the once-strenuous process of personal hygiene.

### Who's using the Support Station?

Children and adults with virtually any amount of weight-bearing ability can use the Support Station. It's ideal wherever clients are too large to be safely transferred to a flat changing surface. And it's proving especially valuable in schools, where it not only speeds hygiene, it also helps kids develop skills for main-stream schooling and independent living.



### How does it make a difference?

- ➊ **Greater dignity and independence:** Putting an end to infant-like, horizontal hygiene can have dramatic emotional benefits. Even people with very limited weight-bearing ability can now participate in their own personal care.
- ➋ **Improved safety for caregivers:** Eliminating transfers to a mat or changing table means less lifting and wrestling with clothing, plus clothing and skin are more accessible – so the process is healthier for everyone.
- ➌ **Practice in strength and skillbuilding:** Using their own strength and skills to transfer and support themselves, people improve weight-bearing and motor skills over time.
- ➍ **Less time out:** Faster than horizontal diapering, the Support Station lets everyone spend less time on hygiene and more time in class, in therapy, or simply enjoying life.

### A user-friendly design

- Compact and easy to install
- Fully adjustable for diverse users of all ages
- Wall-mounted but floor-supported (250 lb max. load)
- Slots, elbow curves and handholds help users pull up
- Padded support strap for extra security
- Padded knee board for leverage during transfer
- Anti-microbial and easy-to-clean surfaces

