



### DIMENSIONS

Width: 13.25" to 20"  
 Height: 3.25"  
 Weight: Approx. 3.0 lbs.

### WEIGHT LIMIT

No weight limit, yet the cushion must be properly sized to the individual.

### MEDICARE CODES

E2622  
 E2623

### LIMITED WARRANTY

Cushion: 24 Months  
 Cover: 6 Months

■ SKIN PROTECTION   □ POSITIONING   ■ ADJUSTABILITY   ■ ADAPTABILITY   ■ COMFORT

# ROHO® MID PROFILE™ Single Compartment Cushion

The MID PROFILE Single Compartment Cushion is made of 3.25" tall interconnected air cells that provide an optimal seating environment to facilitate wound healing. Made with ROHO's DRY FLOATION® Technology, the cushion is intended to conform to an individual's seated shape, minimizing tissue deformation and maintaining blood flow. The MID PROFILE™ Cushion provides a greater contact area with the opportunity for moderate immersion and envelopment to manage the individual's load.

Based on clinical, scientific and/or engineering evidence, this product may be suitable for individuals who:

- Have a history of or are at risk for skin/soft tissue breakdown and/or any category/stage pressure ulcer (including deep tissue injury), skin grafts, or myocutaneous flap surgery
- Currently have skin/soft tissue breakdown and/or any stage pressure ulcer (including deep tissue injury)
- Have normal or impaired sensation
- Require accommodation of pelvic asymmetry, up to approximately 2.25 inches