

ROHO® HIGH PROFILE® Single Compartment Cushion

The HIGH PROFILE Single Compartment Cushion is made of 4.25" tall interconnected air cells that provide an adjustable, low-deformation environment that assists in the healing and prevention of pressure ulcers. The HIGH PROFILE Cushion provides a greater contact area with the opportunity for maximum immersion and envelopment to manage the individual's load.

Based on clinical, scientific and/or engineering evidence, this product may be suitable for individuals who:

- Have a history of or are at risk of skin/soft tissue breakdown and/or any category/stage pressure ulcer (including deep tissue injury), skin grafts, or myocutaneous flap surgery
- Currently have skin/soft tissue breakdown and/or any stage pressure ulcer (including deep tissue injury)
- Have normal or impaired sensation
- Require accommodation of pelvic asymmetry, up to approximately 3 inches



DIMENSIONS

Width: 11.5" to 26.75"
Height: 4.25"
Weight: Approx. 3.25 lbs.

MEDICARE CODES

E2622
E2623

WEIGHT LIMIT

No weight limit, yet the cushion must be properly sized to the individual.

LIMITED WARRANTY

Cushion: 24 Months
Cover: 6 Months

■ SKIN PROTECTION □ POSITIONING ■ ADJUSTABILITY ■ ADAPTABILITY ■ COMFORT