



THE ADAPTOR[®] PAD

The ADAPTOR Pad can be used as an interface on virtually any flat or contoured surface to protect soft tissues. The ADAPTOR Pad is not intended as a seating support surface; it is especially useful when it is cut to fit various surface shapes, such as arm rests, trays, head and foot rests or anywhere uncontrolled movement may need protection.

- Individually sealed air cells
- Can be cut to size
- Absorbs shock
- Protects against shear and friction
- Enhances soft tissue protection

DIMENSIONS

Height: 1"
Width: 9.75"
Depth: 13.5"
Weight: 0.5 lbs.

WEIGHT LIMIT

None, if properly sized to the individual.

LIMITED WARRANTY

24 months