


ROHO®

# AGILITY® Back Supports – Standard series

The ROHO AGILITY back is designed to fit you and your unique needs. It is one of our strongest and easiest to adjust wheelchair back supports. Height, depth, and angle adjustments can be done even when you are seated to reduce the need for transfers. AGILITY backs are quick to install and adaptable to virtually any wheelchair. You can count on them to stay secure and in position to give you the support you need.

AGILITY backs are lightweight and durable thanks to the shell’s aircraft-grade aluminum, and feature a DRY FLOATATION insert to provide you with added stability and comfort. The air insert is fully adjustable and can be tailored to your comfort level and postural needs.



All AGILITY backs can be flipped to accommodate varying body shapes such as wider torsos or broader shoulders.



Lateral air inserts only available on the AGILITY Max Contour.

ROHO

# AGILITY Max

E2615/E2620

The AGILITY Max Contour back support delivers 7” of contour depth to provide significant posterior and lateral trunk support to give you optimal posture and stability. It is the only deep contour back with independent, adjustable, lateral air pads for better positioning and increased comfort.

Dimensions\*    Widths: 13”-20” (33 cm-51 cm)  
                         Depth: With air insert 7” (18 cm)  
                         Lengths: 16”-20” (40.5 cm-51 cm)

*\*For part numbers and full size listings, refer to the order form or visit permobil.com*

ROHO

# AGILITY Minimum

E2611/E2613/E2615/E2620

The AGILITY Minimum Contour back support gives you basic postural support that is lightweight and streamlined to fit your active lifestyle. This back is ideal for individuals who have good trunk control and need adjustability to enhance their range of motion.

Dimensions\*    Widths: 13”-20” (33 cm-51 cm)  
                         Depth: Foam: 2.5” (6.5 cm)  
   Foam/Air: 3” (7.5 cm)  
                         Lengths: 10”-18” (25.5 cm-45.5 cm)

*\*For part numbers and full size listings, refer to the order form or visit permobil.com*



**Zytel® hardware**  
The AGILITY quick release and fixed hardware are made from Zytel – a DuPont™ product which is known for its excellent strength-to-weight ratio, high mechanical strength, and abrasion resistance.

ROHO

# AGILITY Mid

E2615/E2620

The AGILITY Mid Contour back support has 4” of contour depth to provide you with posterior trunk support and moderate lateral support. This back is ideal for individuals who need additional trunk support and may benefit from having an adjustable back, or someone who may need or benefit from increased range of motion.

Dimensions\*    Widths: 13”-20” (33 cm-51 cm)  
                         Depth: With air insert 4” (10 cm)  
                         Lengths: 10”-20” (25.5 cm-51 cm)

*\*For part numbers and full size listings, refer to the order form or visit permobil.com*



**Easy to adjust**  
With the AGILITY back, adjustments can be made up to 4” in height and up to 1” in depth. You can also get up to 15° of anterior or posterior tilt.

**Comfortable**  
AGILITY backs integrate ROHO DRY FLOATATION® technology, an air insert that is fully adjustable and can be tailored to your comfort level and postural needs.