

# ROHO® AirLITE® Cushion

The AirLITE Cushion is designed to improve transfers by having a lower front and middle section, offering positioning capabilities with leg troughs. The lower height of the cushion also contributes to a greater feeling of security in the wheelchair.

A non-adjustable ROHO AIR FLOATATION™ air insert is built into the cushion under the ischials, to help prevent the individual from bottoming out. Softer, more comfortable foam allows the cushion to conform better to the individual's physiology.

Based on clinical, scientific and/or engineering evidence, this product may be suitable for individuals who:

- Are at risk for skin/soft tissue breakdown
- Have normal sensation
- Require increased stability
- Require lower extremity alignment
- Need a lightweight, non-adjustable seating system



\*Cover is required for use.

## DIMENSIONS

Width: 13" to 20"  
Height: Approx. 2" to 3.75"  
Weight: 2.25 lbs.

## MEDICARE CODES

E2605

## WEIGHT LIMIT

300 lbs.

## LIMITED WARRANTY

Cushion: 24 Months  
Cover: 6 Months

SKIN PROTECTION  POSITIONING  ADJUSTABILITY  ADAPTABILITY  COMFORT