



# PRODIGY® Mattress Overlay System

The ROHO® PRODIGY® Mattress Overlay System provides an individualized fit in order to protect skin and soft tissue from breakdown.

A non-powered, air-filled, adjustable, zoned, reactive, full support surface, PRODIGY is intended to conform to an individual's shape to provide skin/soft tissue and deep tissue injury protection. Flexible, interconnected air cells can be adjusted to match the individual's body contour. Each of the three sections can be adjusted separately to address specific, changing needs.

- Zoned - Three sections each 36 in. × 27.5 in.
- Segmented for use in positions other than 180°.
- Hydrostatic distribution of load across the entire contact area.
- Provides 3 in. of immersion.

#### DIMENSIONS

Width: 36"  
Height: 3.5"  
Depth: 81.75"  
Weight: 9.5 lbs.

#### WEIGHT LIMIT

300 lbs. / 136 kg

#### LIMITED WARRANTY

12 Month Limited  
Warranty on soft goods.