



## Hybrid Series

These innovative cushions combine the skin and soft tissue protection from the Shape Fitting Technology™ of an air-filled ROHO DRY FLOATATION™ Insert with the stability of a contoured foam base.

The contoured foam base is designed to allow for increased stability for transfers and provides exceptional comfort, positioning and stability.

- Combination of air and foam for stability and skin/soft tissue protection
- Contoured foam provides lateral and medial upper leg supports to hold your legs in place
- Firm base for transfers
- Supports pressure injury prevention and treatment





## Hybrid Elite SR™ Single Compartment Cushion

This innovative cushion combines the skin and soft tissue protection from the Shape Fitting Technology of an air-filled ROHO DRY FLOATATION® Cushion with stability from a contoured foam base.

The shape-fitting air insert provides critical pressure redistribution around the pelvic area.

---

### IDEAL FOR PEOPLE WHO

Require accommodation of pelvic asymmetry, up to 2 in. / 5 cm

Require pelvic and lower extremity positioning

Require a more stable base for transfers

### DIMENSIONS

Width: 13.75 in. to 23.75 in. / 35 cm to 60 cm

Depth: 14.75 in. to 20.75 in. / 37.5 cm to 53 cm

Height: 4.25 in. / 11 cm

### LIMITED WARRANTY

Cushion: 24 months

Cover: 6 months

### WEIGHT LIMIT

500 lbs. / 225 kg

### MEDICARE CODES

E2622, E2623



## Hybrid Elite Dual Compartment Cushion

Combines the skin and soft tissue protection and Shape Fitting Technology of a ROHO DRY FLOATATION® Cushion with the stability of a contoured foam base.

The dual compartment air insert provides critical pressure redistribution around the pelvis. The two independent, adjustable air compartments let you set inflation levels to reduce or accommodate a side-to-side asymmetry, such as a pelvic obliquity.

---

### IDEAL FOR PEOPLE WHO

Require accommodation or correction of pelvic asymmetry, up to 2 in. / 5 cm

Need side-to-side positioning

Require a more stable environment for transfers

### DIMENSIONS

Width: 14 in. to 20 in. / 35.5 cm to 51 cm

Depth: 14.75 in. to 20.75 in. / 37.5 cm to 53 cm

Height: 4.25 in. / 11 cm

### LIMITED WARRANTY

Cushion: 24 months

Cover: 6 months

### WEIGHT LIMIT

500 lbs. / 225 kg

### MEDICARE CODES

E2622, E2623



## Nexus SPIRIT®

ROHO DRY FLOATATION® Technology combined with the stability of a soft, contoured foam base.

This base provides stability and positioning and support when transferring. The individual air cells envelop your body, providing pressure relief and helping to prevent skin and soft tissue breakdown.

---

### IDEAL FOR PEOPLE WHO

Require accommodation of pelvic asymmetry, up to approx. 1 in. / 2.5 cm

Need transfer stability and support

Need a combination of air and foam for positioning

### DIMENSIONS

Width: 14 in. to 20 in. / 35.5 cm to 51 cm

Depth: 15.5 in. to 19.5 in. / 39.5 cm to 49.5 cm

Air Insert Height: 2.25 in. / 5.5 cm

Total Height: 3.5 in. / 9 cm

### LIMITED WARRANTY

Cushion: 24 months

Cover: 6 months

### WEIGHT LIMIT

No weight limit

### MEDICARE CODES

E2622