



ROHO® Shower / Commode Cushion

The ROHO Shower/Commode Cushion is an adjustable, air-filled, cellular-design, support surface that utilizes DRY FLOATATION® Technology and is intended for short-term cushioning needs on bedside commodes and shower chairs to provide skin and soft tissue protection. The Shower/Commode Cushion is easy to clean and disinfect. Cover not included.

DIMENSIONS

Height: 2.25"
Width: 14.75"
Depth: 14.75"
Weight: 1.5 lbs.

WEIGHT LIMIT

None, if properly sized
to the individual

LIMITED WARRANTY

24 months