


# ROLL-A-RAMP® VAN RAMPS

*The Logical, Affordable Solution to Van Access*



## Why choose a Roll-A-Ramp® Van Ramp?

- Common sense pricing – A fraction of the cost of a converted van or lowered floor conversion
- Keep or purchase a standard factory-model van\* (\*must verify compatibility)
- 6061 Aerospace Aluminum Design - will not rust or corrode
- Ramp stores inside vehicle and is protected from weather elements



**Unique Features**

- All-weather traction
- Folds easy for storage
- 8" to 12" occupied interior space
- Powered and Non-Powered option
- Non-powered can be removed w/ quick pins
- 1,000 lbs load capacity
- Mounts into floor of many standard vans
- \*check for compatibility

## REAR INSTALLATION: ONLY

Transit Connect / ProMaster City / Nissan NV200



## SIDE INSTALLATION: ONLY

Dodge Caravan / Chrysler Town & Country / Toyota Sienna\*

