

Lightweight rollator walker with shopping bag

rollz flex

Walk with confidence every step of the way

- Multi-support options
- Great stability
- Small turning circle
- Soft tires
- Easy to push and steer
- Height adjustable

Designed in the Netherlands



Rollz Flex² keeps you active

At first glance, it does not appear to be a walking aid. This was exactly the intention of the designers: to create an enjoyable design for going outdoors. This rollator walker maneuvers easily and you can walk stable and comfortable, while the bag carries your groceries.

Flexible

Easy adjustable handlebar for different heights.

Practical

The large water-resistant bag can carry up to 44 lbs/6.6 gal. Includes removable inner bag.

Maneuverable

Light steering and easy motion thanks to the small turning cycle.



Convenient

Rest on the leather look seat while having the handlebar as a supportive backrest.

Easy aid

Go over doorsteps and curbs easily with the curb assistant, even with a full loaded bag.

Comfortable

Rollz Flex² allows a comfortable ride with less vibrations thanks to the softer foam-filled tires.

Customize the lightweight rollator walker



Cane holder



Umbrella



Tray



Shopping hooks



Light



Back support

Frame colors



Matt black



Pebble white



Dark purple



Jungle green

Specifications

	medium	small
Rollator weight	17 lb.	15 lb.
User height	5'1" - 6'4"	4'6" - 5'25"
Seat height	24.4"	20.5"
Weight Capacity	275 lb.	275 lb.
Overall width x length	24.5 x 23.6"	24.5 x 23.6"
Width folded	13"	13"
Adjustable height	31.5 - 39.4"	27.5 - 31.9"