

2-in-1 rollator walker and wheelchair  
with air tires

rollz motion  
performance

## Experience comfortable walks even on rough terrain

- Transforms within seconds
- Extra comfortable on rough terrain
- Shock-free and low-vibration
- Easy to push and steer
- Adjustable hand grips and footrests

Designed in the Netherlands



# Have a safe ride with Rollz Motion Performance

The all-terrain rollator walker is very pleasant to use in off road environments such as country paths, over cobblestones or other types of difficult pavement. It provides advanced shock absorption for those who need additional cushioning due to sensitive hands, wrists, shoulders or back problems.

## Soft padding

Armrests with extra thick foam padding to increase comfort.

## Firm seat

Wedge-shaped cushion for a better sitting position.

## Air filled tires

Sturdy air tires enable more comfortable rides on bumpy roads.



## Ergonomic grips

Ergonomically shaped soft handles for optimal grip.

## Visible

Visible in the dark thanks to multiple reflective parts.

## Stable

Walk tall and steadily, while safely supported by a sturdy design.

# Accessorize the all-terrain rollator walker and wheelchair



## Frame color



Jungle green

## Specifications

Rollator weight	25.1 lb.
Wheelchair package weight	7.5 lb.
Maximum weight capacity	275 lb.
User height	5.25 - 6.25 ft
Seat height	21.5"
Handle height range	34 - 39"
Overall width x length	26 x 24.8"
Width folded	11.5"