Parkinson's rollator walker with customizable cues



Regain your walking rhythm

- Laser, sound and vibration cues
- Great stability and comfort
- Advanced and smart cue control
- Easy to push and steer
- Supports various rehabilitation programs





Step out of a freeze with the Rollz Motion Rhythm

This stable rollator walker is developed for people with a disrupted gait pattern or who suffer from freezing episodes. It encourages the start of a rhythmic gait pattern with the help of three cues.

Innovation

Specialized version of the Rollz Motion that addresses people with Parkinson's or a disrupted gait pattern.

Three cues

Laser line on the ground, metronome sound, and vibration in the handlebars.









Laser

Vibration

Sound

Adjustable

Settings can be customized in accordance with the rehabilitation program.

Start the rhythm

The desired cues combination and tempo are easy to set and activate.

Stable

Offers unmatched stability to prevent falls that could occur due to balance problems.

Smartphone app

Allows to create random patterns for preventing habituation.



Accessorize the Parkinson's rollator walker



Frame color



Specifications

Rollator weight	25.55 lb.
Maximum weight capacity	275 lb.
User height	5.25 ft - 6.25 ft
Seat height	21.5"
Handle height range	34 - 39"
Overall width x length	26" x 24.8"
Width folded	12.2"