

Portable ceiling lift comparison

Savaria Monarch


The Monarch portable ceiling lift is one of the lightest on the market weighing only 11 lb, but still lifts the industry standard of 440 lb.

The unique spreader bar includes our industry-leading secure sling mechanism designed to ensure loop attachments cannot be unintentionally removed.

- Ultra-lightweight unit
- Express recharge
- Increased battery capacity
- Unique sling retention spreader bar



Savaria Monarch
portable ceiling lift



Voyager

Maxi Sky 440

Prism P-440

RisePorto 4500

	Savaria Monarch	Voyager	Maxi Sky 440	Prism P-440	RisePorto 4500
Product weight	11 lb	13.3 lb	15.5 lb	11 lb	13.3 lb
Weight capacity	286 lb / 440 lb	440 lb	440 lb	440 lb	308 lb / 452 lb
Speed (in/s) with 264.5 lb	2	1.6	1.6		1
Speed (in/s) with 440 lb	1.8	1.2	1.2		0.7
Travel (in)	90.5	86.6	86.6	78.7	66.9
Soft start/stop	Yes	Yes	Yes	Yes	No
Number of lifts with an avg. charge at 198 lb load	80	26	27	35	50
Lifting clearance	13.5 in	13.5 in	15.8 in	18 in	18 in