

seca 761 & sauna 761 & 760 & 754 & 750

Mechanical personal scales



seca 761



seca sauna 761

- Robust construction.
- Large, easy-to-read dials.
- Low profile and wide platform.
- Precise and durable.



seca 760



seca 754



seca 750

seca 761 & sauna 761 & 760 & 754 & 750:

Attractive precision in a very confined space.

The compact, mechanical flat scales of seca with their durability, functionality and size permit medically precise weighing results in even very limited spaces. The large range gives you the possibility of selecting the right scale to meet your needs in the medical practice and to match your individual style.

seca 761 – the original big dial.

The large, easy-to-read dial gives a truly accurate weight reading. The robust, all-metal body is enhanced by a heavy-duty lever system. Now with a low, wider non-slip rubber platform.

seca sauna 761 – big dial with corrosion protection.

The frame and mechanism of the seca sauna 761 has been optimized for use where there is high humidity – such as, for example, the sauna or swimming baths. The casing has been especially sealed and the mechanics provided with an elaborate corrosion protection. The non-slip platform is made of stainless steel.

seca 760 – big dial in five different colors.

How about a burgundy red, midnight blue, ecru or with white or black housing? Each color comes with its own matching grooved rubber platform mat. The same robust construction and accuracy you expect from seca big dials. The precision movement of this scale with circular display guarantees exact results and the design with a lot of colour and chrome gives the scale a fresh, timeless look.



The large dial makes it easy to read results.

seca 754 – mechanical floor scale with nickelplated casing.

Large, easy-to-read white numbers on black dial. Nickelplated casing. This scale provides the same great precision that big dials are known for. The rugged good looks of the seca 754 give it special appeal in athletic environments.

seca 750 – mechanical floor scale.

The mechanical circular dial scale seca 750 possesses a rugged durability with a big dial, at a lower price thanks to 21st century design. A lower profile and wider platform enhance utility. Available in lbs or lbs/kg versions.



seca 761 & sauna 761 & 760 & 754 & 750

Technical data

seca 761:

- Capacity: 330 lbs / 150 kg
- Graduation: 1 lbs / 1,000 g
- Dimensions (WxHxD):
11.9 x 4.6 x 18.5" /
303 x 118 x 470 mm
- Platform (WxD): 11 x 11" /
279.4 x 279.4 mm
- Weight: 7.7 lbs / 3.5 kg

seca sauna 761:

- Capacity: 330 lbs / 150 kg
- Graduation: 1 lbs / 1,000 g
- Dimensions (WxHxD):
11.9 x 4.6 x 18.5" /
303 x 118 x 470 mm
- Platform (WxD): 11 x 11" /
279.4 x 279.4 mm
- Weight: 7.7 lbs / 3.5 kg

seca 760:

- Capacity: 330 lbs / 150 kg
- Graduation: 1 lbs / 1,000 g
- Dimensions (WxHxD):
11.9 x 4.6 x 18.5" /
303 x 118 x 470 mm
- Platform (WxD): 11 x 11" /
279.4 x 279.4 mm
- Weight: 7.7 lbs / 3.5 kg

seca 754:

- Capacity: 330 lbs / 150 kg
- Graduation: 1 lbs / 1,000 g
- Dimensions (WxHxD):
12.5 x 4.6 x 18.8" /
317 x 117 x 477 mm
- Platform: (WxD) 11 x 11.5" /
279.4 x 292.1 mm
- Weight: 7.7 lbs / 3.5 kg

seca 750:

- Capacity: 330 lbs / 150 kg
- Graduation: 1 lbs / 1,000 g
- Dimensions (WxHxD):
12.5 x 4.6 x 18.8" /
317 x 117 x 477 mm
- Platform (WxD): 11 x 11.5" /
279.4 x 292.1 mm
- Weight: 7.7 lbs / 3.5 kg