


SB6c

Dimensions = Inches




Product features:

- Weight capacity - 300lb
- Rolls over most standard toilets
- Low maintenance
- Arms - lockable / flip-back
- Height adjustable
- Footrests - removable and swing-away, height adjustable
- No rust aluminum frame (Grade 6061 T6)
- Stainless steel parts and hardware (Grade 304)
- Safety belts - lap or chest belts plus calf strap
- Caster Wheels - 4" all locking
- Commode system with deflector splash guard
- Tool-free assembly

Accessory items

- Refer to the accessory pages or on our website



SB6c ROLL-IN BUDDY LITE

For roll in showers



The SB6c mobile shower chair from Showerbuddy is designed for accessible or roll-in showers. It can roll straight over most toilets and into the shower, reducing manual handling and transfers for caregivers.

The SB6c has a 450mm wide foam seat cushion and can also be used as a commode chair. The SB6c shower chair comes with 4 locking high quality caster wheels, a commode system and numerous accessories to choose from.

The SB6c Roll-InBuddy Lite is lightweight, durable and has a lifetime warranty on the frame.

Product Features:

- Weight capacity - 136kg
- Rolls over most standard toilets
- Arms - lockable / flip-back
- Height adjustable
- Footrests - removable and swing-away, height adjustable
- No rust aluminum frame (Grade 6061 T6)
- Stainless steel parts and hardware (Grade 304)
- Safety belts - lap or chest belts plus calf strap
- Caster Wheels -100mm all locking
- Commode system with deflector splash guard
- Tool-free assembly

Product Dimensions:

- Chair width - over caster wheels 505mm
- Chair Width – over arms 565mm
- Chair Width – Inside Arms 470mm
- Chair Length - including footrest 890mm (630mm with footrests removed)
- Backrest Width 448mm
- Backrest Height from Cushion to top 380mm
- Toilet Height Clearance (maximum) 490mm
- Toilet width clearance at narrowest 365mm (at casters)
- Toilet clearance can depend upon the shape of the toilet pan
- Seat Cushion Dimensions 460 W x 430mm D
- Seat Height from Floor 510mm to 585mm (Adjustable)
- Seat Cushion Opening 175mm W x 335mm L
- Box Dimensions L 580mm x W 570mm x H 500mm
- Box Gross Weight 15.8kg

