

# SB6w

## ROLL-INBUDDY SOLO

For roll-in showers



The SB6w mobile shower chair from Showerbuddy is designed for accessible or roll-in showers and for users who can self-propel. It can roll straight over most toilets and into the shower, reducing transfers for the user. The SB6w has a 450mm wide soft foam seat cushion that is height adjustable and can also be used as a commode chair.


The SB6w shower chair comes with 2 x 24" wheels with hand rims and 2 x locking high quality caster wheels, a commode system and numerous accessories to choose from.

The SB6w Roll-InBuddy Solo is lightweight, durable and has a lifetime warranty on the frame.



# SB6w

Dimensions = Inches




## Product features:

- Weight capacity - 300lb
- Rolls over most standard toilets
- Low maintenance
- Arms - Lockable / flip-back
- Height adjustable
- Footrests - removable and swing-away, height adjustable
- No rust aluminum frame (Grade 6061 T6)
- Stainless steel parts and hardware (Grade 304)
- Safety belts - lap or chest belts plus calf strap
- Caster Wheels - 4" all locking
- Commode system with deflector splash guard
- Soft blow-mold cushion
- Tool-free assembly

## Accessory items

- Refer to the accessory pages or on our website



# SB6w ROLL-IN BUDDY SOLO

## For roll in showers

The SB6w mobile shower chair from Showerbuddy is designed for accessible or roll-in showers and for users who can self-propel. It can roll straight over most toilets and into the shower, reducing transfers for the user. The SB6w has a 450mm wide soft foam seat cushion that is height adjustable and can also be used as a commode chair.

The SB6w shower chair comes with 2 x 24" wheels with hand rims and 2 x locking high quality caster wheels, a commode system and numerous accessories to choose from.

The SB6w Roll-InBuddy Solo is lightweight, durable and has a lifetime warranty on the frame.



### Product Features:

- Weight capacity - 136kg
- Rolls over most standard toilets
- Low maintenance Arms - Lockable / flip-back
- Height adjustable
- Footrests - removable and swing-away, height adjustable
- No rust aluminum frame (Grade 6061 T6)
- Stainless steel parts and hardware (Grade 304)
- Safety belts - lap or chest belts plus calf strap
- Caster Wheels - 100mm all locking
- Commode system with deflector splash guard
- Soft blow-mold cushion
- Tool-free assembly



### Product Dimensions:

- Chair Width - over wheelchair wheels 635mm
- Chair Width - over arms 565mm
- Chair Width - Inside Arms 470mm
- Chair Length - including footrest 1060mm (820mm with footrests removed)
- Backrest Width 448mm
- Backrest Height from Cushion to top 380mm
- Toilet Height Clearance (maximum) 495mm
- Toilet width clearance at narrowest 365mm (at casters)
- Toilet clearance can depend upon the shape of the toilet pan
- Seat Cushion Dimensions 460 W x 430mm D
- Seat Height from Floor 570mm to 595mm (Adjustable)
- Seat Cushion Opening 175mm W x 335mm L
- Box Dimensions L 660mm x W 640mm x H 540mm
- Box Gross Weight 20.8kg

