

pediatric positioning

Skillbuilders® strap wedge



- two 5" wide positioning straps provide support when positioning the child
- without the removable straps, the deluxe wedge functions just like the standard Skillbuilders® wedge
- used with optional leg abductor wedge (30-1018 shown on page 219)

	L x W x H		incline
30-1020	22" x 20" x 4"	small	11°
30-1021	22" x 20" x 6"	small	16°
30-1022	22" x 20" x 8"	small	20°
30-1023	22" x 20" x 10"	small	25°
30-1024	26" x 24" x 6"	medium	9°
30-1025	26" x 24" x 8"	medium	13°
30-1026	26" x 24" x 10"	medium	17°
30-1027	26" x 24" x 12"	medium	22°

Skillbuilders® adjustable wedge system



- arms and hands can move freely while lying on this adjustable wedge
- both the front and rear legs can be adjusted to the proper height with single button locks
- adjustable: front 2-15", rear 2-7"
- removable abductor adjusts 4-8"
- side rail for carrying, securing positioning straps
- wedge folds flat
- vinyl covered foam padded top

31-1050 adjustable wedge (22"W x 36"L)

Skillbuilders® bi-height wedge



- use at higher or lower height (side 1 or side 2)
- high molded sides give lateral support
- positioning strap provides support
- accommodates abductor wedge (30-1018 shown on page 219)

bi-height wedge					
	L x W	side 1	incline	side 2	incline
30-1240	22" x 18"	4"	11°	6"	16°
30-1241	22" x 18"	8"	20°	10"	25°
30-1242	26" x 22"	6"	9°	8"	13°
30-1243	26" x 22"	10"	17°	12"	22°

Versa-Form® positioning pillows



- provides customized, semi-permanent (regular: 1 month, plus: 3 months) positioning
- shape/reshape in seconds by opening air valve and molding to desired shape and then evacuating air with vacuum pump

flat size	Versa-Form®	Versa-Form® Plus	covers
16 x 20"	30-2024	30-2004	30-2034
22 x 12"	30-2020	30-2000	30-2030
22 x 26"	30-2021	30-2001	30-2031
22 x 34"	30-2022	30-2002	30-2032
32 x 40"	30-2023	30-2003	30-2033

vacuum pump

- single handed operation allows the user to form the pillow

30-2040 pump

