


Skillbuilders® swing seats



seat	inner dimensions	capacity
small	20"L x 12"W x 9"D	75 lb
medium	23"L x 15"W x 10"D	150 lb
large	31"L x 20"W x 14"D	200 lb

- plastic seat is available with and without a built-in abductor (pommel)
- safety harness and crotch strap standard for proper positioning
- adjustable suspension mechanism will fit any frame up to ten feet
- optional headrest, legrest, seat cover and swing tie-down straps available

- complete full support swing seat comes with pommel, headrest and legrest, harness and crotch strap
- small is for child, medium for teenager and large for adult

	seat with pommel	
	with rope	with chain
small	30-1630	30-1590
medium	30-1631	30-1591
large	30-1632	30-1592

	seat without pommel	
	with rope	with chain
small	30-1634	30-1593
medium	30-1635	30-1594
large	30-1636	30-1595

	complete swing seat	
	with rope	with chain
small	30-1640	30-1596
medium	30-1641	30-1597
large	30-1642	30-1598


bouncer accessory



- bounce adapter adds a new dimension of vestibular challenge to swinging
- child can bounce while swinging back-and-forth or use without swinging
- 125 lb working load each bouncer

30-1692	bouncer, each
30-1692-2	bouncer, pair

rotator accessories



- rotator adapter adds a new dimension and vestibular challenge to swinging
- | | |
|---------|----------------------|
| 30-1695 | single attachment |
| 30-1697 | quadruple attachment |

swing seat accessories



headrest	
30-1633	small/medium
30-1629	large

legrest	
30-1643	small/medium
30-1644	large

covers for seat with pommel (not shown)	
30-1645	small
30-1646	medium
30-1647	large

replacement seat harnesses (not shown)	
30-1630H	small
30-1631H	medium
30-1632H	large

- **tie-down retention straps (not shown)**
 - anchor swing to ground to limit motion
- | | |
|---------|-----------------|
| 30-1637 | straps (2 each) |
|---------|-----------------|