

# swing frames for vestibular training and play therapy

## Vestibular A-frame complete set

- durable steel free-standing frame
- 12'10"L x 6'6"W x 7'6"H
- 500 lb capacity; 300 lb working load
- set includes: frame, 4-point spinner, log swing, circle platform swing, flying saucer swing, therapy net swing with board, vertical bouncer, 2 height adjustment systems, 2 safety rotational devices, 4 safety snaps and 10' x 5' rainbow mat

30-1600 A-frame complete set



## swing frames

### indoor/outdoor frame

- portable indoor/outdoor frame made with galvanized-steel pipe construction
- 83"L x 70"W x 77"H
- 250 lb capacity

A 30-1671 portable frame

### indoor frame

- with green powder-coat finish
- 96"L x 96"W x 96"H
- 500 lb capacity

B 30-1672 indoor frame

