



# SleepSafe® Bed

Quote Request  Order Request

(Please Check One)

Company Name

Requested By

Ship To Address

Phone

Fax

City

State

Zip

Email

Billing Address

Your Client Name for Reference

City

State

Zip

Your Purchase Order # (Required at Order)

SSB Quote #

(Please CHOOSE ONLY ONE choice per numbered category – EXCEPT for Optional Accessories and BASIC Bed ordering)

<b>1. Safety Rail Height</b>	SS	<input type="checkbox"/>	LOW - SleepSafe® Bed
	S2	<input type="checkbox"/>	MEDIUM - SleepSafe® II Bed
	SR	<input type="checkbox"/>	TALL - SleepSafer® Bed
	OEXT	<input type="checkbox"/>	TALL - SleepSafer® / OPEN Extension with Flip-Up Rail    Add optional padding to extension
	CEXT	<input type="checkbox"/>	TALL - SleepSafer® / CLOSED Extension with Flip-Up Rail    Add optional padding to extension

<b>2. Foundation</b>	BB or BX	<input type="checkbox"/>	FIXED No movement.    Bunkie Board (5" more safety rail height) <u>OR</u> Box Spring
	MP	<input type="checkbox"/>	MANUAL ARTICULATION    Elevates head and knee only, with hand cranks
	ME	<input type="checkbox"/>	MANUAL ELECTRIC COMBO    Elevates mattress height with hand crank <i>plus</i> head & knee with electric remote
	MF	<input type="checkbox"/>	MANUAL HI-LO    Elevates mattress height only, with hand crank
	MH	<input type="checkbox"/>	MANUAL HI-LO COMBO    Elevates mattress height with hand crank <i>plus</i> head & knee with hand cranks
	PL	<input type="checkbox"/>	ELECTRIC ARTICULATION    Elevates head and knee only, with electric remote
	HL	<input type="checkbox"/>	ELECTRIC HI-LO COMBO    Elevates mattress height plus head & knee with electric remote

<b>3. Type of Finish</b>	OA	<input type="checkbox"/>	OAK (Natural Oak with clear finish; highlights natural wood grain. Cream trim)
	MA	<input type="checkbox"/>	MAPLE (Natural Maple with clear finish; very light with little wood grain. Cream trim)
	AL	<input type="checkbox"/>	ALDER (Rustic birch with random knots and straight grain. Moderately light. Clear finish. Gray trim)
	CH	<input type="checkbox"/>	CHERRY (Natural Cherry with clear finish; NOT dark red stain. Gray trim)
	MG	<input type="checkbox"/>	MAHOGANY (Natural Mahogany with a clear finish. Will darken with age. Gray trim)
	MC	<input type="checkbox"/>	SOLID COLOR    White    Black    Blue    Green    Red    Yellow    Orange    Purple    Pink
	MC	<input type="checkbox"/>	MULTI COLOR    For multi-color option, please attach color selection form. (All have maple posts and white trim)

<b>4. Bed Size</b>	T	<input type="checkbox"/>	TWIN - (Mattress size: 38" x 77" without pads; 36" x 77" with pads)
	F	<input type="checkbox"/>	FULL - (Mattress size: 51" x 77" without pads; 50" x 77" with pads)
	Q	<input type="checkbox"/>	QUEEN - (Mattress size: 61" x 77" without pads; 60" x 77" with pads)

<b>5. Optional Accessories</b> (Select all that apply)	PD	<input type="checkbox"/>	PADDING    Over Windows    Around Windows    Cream    Gray    Light Blue    Taupe
	WNHB	<input type="checkbox"/>	HEADBOARD WINDOW    Solid Clear    Screen Mesh
	WNFB	<input type="checkbox"/>	FOOTBOARD WINDOW    Solid Clear    Screen Mesh
	IVPL	<input type="checkbox"/>	IV POLE (Comes as a Kit so that it can be attached anywhere on bed)
	TCCO	<input type="checkbox"/>	TUBING CHANNEL ACCESS CUT OUT (Specify    L    R    C in Headboard or    L    R    C in Footboard)
	DV	<input type="checkbox"/>	DUAL VIEW (2 additional safety rails; allows access from both sides; Tall Bed only)
	SFVT	<input type="checkbox"/>	SCREEN MESH WINDOWS (Order to replace solid clear windows on safety rails)

<b>BASIC TWIN Bed Ordering</b> <i>Use this area <u>only</u> to order this bed.</i>	<input type="checkbox"/>	FIXED BUNKIE BOARD    ELECTRIC ARTICULATION    MANUAL ARTICULATION /    Oak    Maple    Alder
	<input type="checkbox"/>	PADDING    Over Windows    Around Windows    Cream    Gray    Light Blue    Taupe