

# SleepSafe® Bed Model Measurements

Use this chart to compare the different measurement specifications from model to model.

<i>SleepSafe® Bed Specifications</i>	<i>Low Bed SleepSafe®</i>	<i>Medium Bed SleepSafe® II</i>	<i>Tall Bed / Extension Bed SleepSafer®</i>
<b>FOUNDATION SAFETY RAIL HEIGHTS</b> <i>The distance from the top of mattress to the top of the safety rail.</i>			
<b>Fixed w/Box Spring</b> (Frame attachment choice)	12" or 8¾"	17¾", 14" or 10"	31" or 24¾" / 47¾" or 41½"
<b>Fixed w/Bunkie Board</b> (Frame attachment choice)	17¼" or 14"	27" or 21¾"	36¼" / 53"
<b>Manual Articulation</b> (Frame attachment choice)	15" or 11¾"	21", 17" or 13¼"	34" or 27¾" / 50¾" or 44½"
<b>Manual Electric Combo</b> (LO to HI Settings based on frame attachment)	LO 11½" to EVEN with rail	LO 21¼" to HI 9½" or LO 16" to HI 4¼"	LO 30½" to HI 18¾" / LO 47½" to HI 35¾"
<b>Manual Hi-Lo</b> (LO to HI Settings based on frame attachment)	LO 11½" to EVEN with rail	LO 21¼" to HI 9½" or LO 16" to HI 4¼"	LO 30½" to HI 18¾" / LO 47½" to HI 35¾"
<b>Manual Hi-Lo Combo</b> (LO to HI Settings based on frame attachment)	LO 11½" to EVEN with rail	LO 21¼" to HI 9½" or LO 16" to HI 4¼"	LO 30½" to HI 18¾" / LO 47½" to HI 35¾"
<b>Electric Articulation</b> (Frame attachment choice)	15" or 11¾"	21", 17" or 13¼"	34" or 27¾" / 50¾" or 44½"
<b>Electric Hi-Lo Combo</b> (LO to HI Settings based on frame attachment)	LO 11½" to HI 1¼" below rail	LO 21¼" to HI 11" or LO 16" to HI 5¾"	LO 30½" to HI 20¾" / LO 47½" to HI 37¼"
<b>FOUNDATION MATTRESS TO FLOOR HEIGHTS</b> <i>The distance from the top of the mattress to the floor.</i>			
<b>Fixed w/Box Spring</b> (Frame attachment choice)	23¾" or 27"	28", 32" or 35½"	23¾" or 30"
<b>Fixed w/Bunkie Board</b> (Frame attachment choice)	18½" or 21¾"	18¾" or 24"	18½"
<b>Manual Articulation</b> (Frame attachment choice)	20¾" or 24"	25", 28¾" or 32½"	20¾" or 27"
<b>Manual Electric Combo</b> (LO to HI Settings based on frame attachment)	LO 24¼" to HI 36"	LO 24½" to HI 36¼" or LO 29¾" to HI 41½"	LO 24¼" to HI 36"
<b>Manual Hi-Lo</b> (LO to HI Settings based on frame attachment)	LO 24¼" to HI 36"	LO 24½" to HI 36¼" or LO 29¾" to HI 41½"	LO 24¼" to HI 36"
<b>Manual Hi-Lo Combo</b> (LO to HI Settings based on frame attachment)	LO 24¼" to HI 36"	LO 24½" to HI 36¼" or LO 29¾" to HI 41½"	LO 24¼" to HI 36"
<b>Electric Articulation</b> (Frame attachment choice)	20¾" or 24"	25", 28¾" or 32½"	20¾" or 27"
<b>Electric Hi-Lo Combo</b> (LO to HI Settings based on frame attachment)	LO 24¼" to HI 34½"	LO 24½" to HI 34¾" or LO 29¾" to HI 40"	LO 24¼" to HI 34½"
<b>FLOOR TO RAIL HEIGHT</b> <i>The distance from the floor to the top of the rail.</i>	35¾"	45¾"	54¾" / 72¾"
<b>OVERALL OUTSIDE DIMENSIONS</b>			
<b>Twin</b>	42¼" W x 80½" L x 40¾" H	42¼" W x 80½" L x 50¾" H	42¼" W x 80¼" L x 56" H (72¾" H w/Ext.)
<b>Full</b>	56" W x 80" L x 40¾" H	56" W x 80" L x 50¾" H	56" W x 80" L x 57" H (72¾" H w/Ext.)
<b>Queen</b>	65" W x 80" L x 40¾" H	65" W x 80" L x 50¾" H	65" W x 80" L x 57" H (72¾" H w/Ext.)
<b>MATTRESS SIZE</b> <i>Width without pads / Width with pads / Length</i>			
<b>Twin / Full / Queen</b>	38" W / 36" W / 77" L	51" W / 50" W / 77" L	61" W / 60" W / 77" L
<b>UNDER BED FLOOR CLEARANCE</b>			
<b>Standard Models</b>	8"	8"	7¾"
<b>TRANSFER HEIGHT</b> <i>The distance from top edge of the folded down safety rail to the floor.</i>	22¼" with Safety Rail On 20¾" with Safety Rail Off	32" with Safety Rail On 30½" with Safety Rail Off	26" with Safety Rail On 24½" with Safety Rail Off



Photo Reference Guide for Bed Measurements

