


OUR MOST ECONOMICAL SAFETY BED!

SleepSafe® BASIC Safety Bed



The New SleepSafe® BASIC Safety Bed Features:

1) Low 21¼" Transfer Height

May provide easier transfer to wheelchair

2) Foundation Options Include:

- TWIN Size - Fixed with Bunkie Board
- TWIN Size - Electric Head & Knee Articulation
- TWIN Size - Manual Head & Knee Articulation

4) Alder, Maple and Oak Finishes Only

5) Padding Option is Available

Other Specifications:

SAFETY RAIL HEIGHT:

- With Fixed Bunkie Board: 22" or 19"
- With Electric Articulating Foundation: 20" or 16½"

MATTRESS TO FLOOR HEIGHT:

- With Fixed Bunkie Board: 17¾" or 21"
- With Electric Articulating Foundation: 23¼" or 20"



SleepSafe® - BASIC Twin Bed in Alder with Safety Rail in the down position.

SleepSafe® Beds is a domestic USA manufacturer of custom safety beds, featuring removable safety side rails, designed to address the issues of falls and entrapment for those with special needs. Models include SleepSafe®-Low Bed, SleepSafe® II- Medium Bed, and SleepSafer®- Tall Bed, offering progressively more safety rail to mattress height.



SleepSafeBed

SleepSafe® Beds are built to order, offering twin, full or queen size; fixed, articulating or HiLo electric or manual foundations; padding; wood or color finishes, and other helpful options. SleepSafe® Beds help those with special needs get a safe, restful sleep with peace of mind for their caregivers.

FDA Registered