

# hardrock

Colour range



## Robust and durable

The Hardrock is an extremely durable and robust chair that has been specifically designed for children who are prone to rock repetitively. The Hardrock provides moderate levels of postural support and is suitable for both school and home environments.

## Tailoring to the child's needs

It can be a challenge to engage children with complex learning difficulties in educational and everyday activities. The Hardrock tackles this issue and offers a positive solution. It provides a safe environment and one that is conducive to learning.

## Posture matters

The Hardrock uses an anti-thrust system and anterior pad to provide secure postural management within clearly defined boundaries. Thoracic supports and stabilisers are available for children who exhibit sideways motion.

## Hardrock technical data

all dimensions in mm unless otherwise stated

Code	Size	Seat height	Width	Seat depth	Weight limit (kg)	Back height
8202	2	270-410	350	300-415	60	390
8203	3	290-430	430	370-485	70	470
8204	4	375-525	470	400-550	90	575

# hardrock

## accessories

Sizes 2 & 3 of the Hardrock chair come complete with adjustable armrests, anti-thrust seat cushion, tray, buckle pelvic strap with pad and transporter base (consists of skis with transporter wheels and integrated footboard).

Size 4 of the Hardrock chair comes complete with adjustable armrests, anti-thrust seat cushion, tray, a buckle pelvic strap with pad and mobile base with push handle.



### Mobile base with push handle

CODE MB010 - MB012

Enables child to be easily moved whilst seated

Comes as standard for size 4 only



### Footplate & heel bar

CODE FP017 - FP019

Provides extra foot support

Straps and sandals can be added



### Stabilisers<sup>1</sup>

CODE SB001 - SB002

Prevents sideways rocking

### Extended side gables<sup>1</sup>

CODE EG001 - EG002

Keeps child's legs within the chair



### Pommel<sup>1</sup>

CODE PL020 - PL021

Encourages hip abduction



### 4-point pelvic strap<sup>1</sup>

CODE HN074 - HN076

Firm support to stabilise hip position



### Increased height backrest (extra 100mm)<sup>2</sup>

CODE BK155 - BK157

Provides extra back support for taller children

<sup>1</sup> Must be factory fitted at time the chair is ordered

<sup>2</sup> Available with Slatted and Plain backrests

# hardrock

## accessories

Colour range



### Tray handrail full width

CODE GB010 - GB012

Helps child to stabilise upper limb  
Additionally can be used to attach toys



### Padded tray insert (fitted top)<sup>1</sup>

CODE TR035 - TR037

Provides comfort and protection



### Padded tray insert<sup>1</sup>

CODE TR031 - TR033

Provides extra comfort and protection  
Fitted underneath the tray



### Front tray retaining strap<sup>1</sup>

CODE ST053

Secures the tray and keeps the tray  
in the desired position



### Rear tray retaining strap<sup>1</sup>

CODE ST054

Secures the tray and keeps the tray  
in the desired position



### Side tray retaining strap<sup>1</sup>

CODE ST055

Secures the tray and keeps the tray  
in the desired position



### Moulded pelvic support<sup>2</sup>

CODE SP017 - SP019

Provides snug contoured support



### Trunk supports with straps<sup>2</sup>

CODE SP010 - SP011

Provides additional trunk support



### Samba style headrest

CODE HT014 - HT015

Provides extra support for the head  
Can be shaped by adjusting the sides  
Available in sizes 1 & 2

<sup>1</sup> Must be factory fitted at time the product is ordered

<sup>2</sup> Requires Slatted backrest

# hardrock

Colour range



## Robust and durable

The Hardrock is an extremely durable and robust chair that has been specifically designed for children who are prone to rock repetitively. The Hardrock provides moderate levels of postural support and is suitable for both school and home environments.

## Tailoring to the child's needs

It can be a challenge to engage children with complex learning difficulties in educational and everyday activities. The Hardrock tackles this issue and offers a positive solution. It provides a safe environment and one that is conducive to learning.

## Posture matters

The Hardrock uses an anti-thrust system and anterior pad to provide secure postural management within clearly defined boundaries. Thoracic supports and stabilisers are available for children who exhibit sideways motion.

## Hardrock technical data

all dimensions in mm unless otherwise stated

Code	Size	Seat height	Width	Seat depth	Weight limit (kg)	Back height
8202	2	270-410	350	300-415	60	390
8203	3	290-430	430	370-485	70	470
8204	4	375-525	470	400-550	90	575

# hardrock

## accessories

Sizes 2 & 3 of the Hardrock chair come complete with adjustable armrests, anti-thrust seat cushion, tray, buckle pelvic strap with pad and transporter base (consists of skis with transporter wheels and integrated footboard).

Size 4 of the Hardrock chair comes complete with adjustable armrests, anti-thrust seat cushion, tray, a buckle pelvic strap with pad and mobile base with push handle.



### Mobile base with push handle

CODE MB010 - MB012

Enables child to be easily moved whilst seated

Comes as standard for size 4 only



### Footplate & heel bar

CODE FP017 - FP019

Provides extra foot support

Straps and sandals can be added



### Stabilisers<sup>1</sup>

CODE SB001 - SB002

Prevents sideways rocking

### Extended side gables<sup>1</sup>

CODE EG001 - EG002

Keeps child's legs within the chair



### Pommel<sup>1</sup>

CODE PL020 - PL021

Encourages hip abduction



### 4-point pelvic strap<sup>1</sup>

CODE HN074 - HN076

Firm support to stabilise hip position



### Increased height backrest (extra 100mm)<sup>2</sup>

CODE BK155 - BK157

Provides extra back support for taller children

<sup>1</sup> Must be factory fitted at time the chair is ordered

<sup>2</sup> Available with Slatted and Plain backrests

# hardrock

## accessories

Colour range



### Tray handrail full width

CODE GB010 - GB012

Helps child to stabilise upper limb  
Additionally can be used to attach toys



### Padded tray insert (fitted top)<sup>1</sup>

CODE TR035 - TR037

Provides comfort and protection



### Padded tray insert<sup>1</sup>

CODE TR031 - TR033

Provides extra comfort and protection  
Fitted underneath the tray



### Front tray retaining strap<sup>1</sup>

CODE ST053

Secures the tray and keeps the tray  
in the desired position



### Rear tray retaining strap<sup>1</sup>

CODE ST054

Secures the tray and keeps the tray  
in the desired position



### Side tray retaining strap<sup>1</sup>

CODE ST055

Secures the tray and keeps the tray  
in the desired position



### Moulded pelvic support<sup>2</sup>

CODE SP017 - SP019

Provides snug contoured support



### Trunk supports with straps<sup>2</sup>

CODE SP010 - SP011

Provides additional trunk support



### Samba style headrest

CODE HT014 - HT015

Provides extra support for the head  
Can be shaped by adjusting the sides  
Available in sizes 1 & 2

<sup>1</sup> Must be factory fitted at time the product is ordered

<sup>2</sup> Requires Slatted backrest