

bath chair



The Bath Chair comes with an adjustable harness and adjustable pommel as standard.

Colour range

White

Bath Chair technical data

all dimensions in mm unless otherwise stated

Size	1	2
Code	7607	7609
Seat height	35	35
Seat width	350	350
Seat depth	250	325
Back height	470	560
User weight (kg)	35	35
Shipping weight (kg)	5	8

Slimline and safe

Bath time for children is an opportunity to play with siblings and parents. Our Bath Chair's slimline design gives your child more space and freedom to bathe without restrictions. Handling the chair is also easy due to its lightweight frame.

Adjustable

Our Bath Chair has an easy adjustable pommel so that legs can be positioned slightly apart for stability. The harness also offers various adjustments to suit your child's needs.

bath chair accessories



Telescopic handrail

CODE HR012

Length range: 470-700mm

Shipping weight: 1kg

Allows child to hold and maintain position