

# bath corner chair

**Colour range**
 White
**More support**

Our Bath Corner Chair is designed for children who require a little extra help with sitting balance whilst bathing. This is achieved by its angled back and hand rail, both of which provide additional security.

**Adjustable**

Like the Bath Chair, our Bath Corner Chair provides numerous levels of adjustment. The pommel and pelvic strap can be adjusted to suit each individual child's needs.

**Compact**

Lightweight and easy to clean, the Bath Corner Chair has a perfectly compact design. This allows it to sit low in the bath enabling your child to fully extend their legs in complete comfort. Its size also enables it to fit a wide range of baths.

**Bath Corner Chair technical data**

all dimensions in mm unless otherwise stated

Size	1	2
Code	7511	7512
Seat height	75	75
Seat width	340	340
Seat depth	280	330
Back height	350	460
User weight (kg)	30	35
Shipping weight (kg)	3	4

# bath corner chair



## accessories

### Colour range

White



The Bath Corner Chair comes complete with velcro pelvic strap, adjustable pommel, telescopic back rail and telescopic handrail as standard.



### Telescopic handrail

CODE HR012

(pictured on the Bath Chair)

Holding the handrail helps to maintain seating balance



### Bath back pad

CODE PD025

Provides extra comfort



"The Bath Corner Chair gives my child extra support and helps me to bathe her."

Christine, Mum, Reading