

# bath corner chair



The Bath Corner Chair comes complete with velcro pelvic strap, adjustable pommel, telescopic back rail and telescopic handrail as standard.

## More support

Our Bath Corner Chair is designed for children who require a little extra help with sitting balance when bathing. This is achieved by its angled back and hand rail, both of which provide additional security.

## Adjustable

Like the Bath Chair, our Bath Corner Chair provides numerous levels of adjustment. The pommel and pelvic strap can be adjusted to suit each individual child's needs.

## Compact

Lightweight and easy to clean, the Bath Corner Chair has a perfectly compact design. This allows it to sit low in the bath enabling your child to fully extend their legs in complete comfort. Its size also enables it to fit a wide range of baths.

## Bath Corner Chair technical data

all dimensions in inches unless otherwise stated

Size	1	2
Item #	7511	7512
Seat height	3.0	3.0
Seat width	13.4	13.4
Seat depth	11.0	13.0
Back height	13.8	18.1
4 Pommel settings (distance to back of pommel)	4.3 / 5.1 / 6.3 / 7.5	5.5 / 6.7 / 7.9 / 9.1
Max user weight (lbs)	66.1	77.2
Product weight (lbs)	6.6	8.8

