

# combi



## Progressive toileting

Our Combi toileting chair is a great starter chair. Its design enables it to be used as a supportive chair while at the same time providing numerous features to assist your child with the everyday toileting needs.

## Providing solutions

The Combi toileting chair's design also incorporates a potty for younger children. This can easily be removed to provide a commode seat for children who have progressed their toilet training.

## Adjustability

Good levels of adjustability are offered with the Combi toileting chair so that a wide range of needs can be met. Seat height and depth can be adjusted to suit the needs of your child.

The Combi Chair comes complete with potty, plastic commode seat, upholstered seat and velcro pelvic strap as standard.

## Combi technical data

all dimensions in inches unless otherwise stated

Size	1	2	3	4
Item #	7551	7552	7553	7554
Seat to floor height*	8.5 - 11.4	10.8 - 14.8	11.8 - 15.7	13.4 - 17.7
Seat width	11.4	12.8	14.4	17.1
Seat depth	8.1 - 11.4	9.6 - 14.0	12.2 - 15.2	13.8 - 17.3
Back height	12.8	14.6	15.7	19.7
Max user weight (lbs)	77.2	99.2	121.3	143.3
Product weight (lbs)	24.3	28.7	30.9	33.1

\*Seat height to floor is to the top of upholstered seat (deduct 1in for the height of the Potty and deduct 1.6in for the height of the Commode seat)

# combi accessories



## Skis & full footrest<sup>1</sup>

ITEM # SK037 - SK040

Provides foot support and offers greater stability



## Reversible footboard<sup>1</sup>

ITEM # FB023 - FB026

Provides foot support  
Prevents unwanted chair movement



## Independent adjustable footrest<sup>1</sup>

ITEM # 4187

Portable, height and angle adjustable  
Solid support for the feet



## Skis

ITEM # SK051 - SK054

Provides greater stability



## Footrest & straps<sup>1</sup>

ITEM # FP002 - FP005

Ensures feet are firmly in position



## Splash guard

ITEM # SG001 - SG002

Assists in leg abduction and keeps toileting contained  
Requires Toggle handrail



## Pelvic cushions single or double thickness (pair)<sup>2</sup>

Provides hip and thigh support  
Available in various sizes



## Toggle handrail

ITEM # HR014 - HR017

Allows child to hold and maintain balance and to sit with minimal support



## Tray

ITEM # TR003 - TR006

Provides trunk support  
Additionally can be used for activities such as playing

<sup>1</sup> Only available in conjunction with a Mobile dolly base, Reversible footboard or Skis with half & full footrest

<sup>2</sup> Available as double cushions (pair) size 1+1in, 1.5+1in & 2+1in and singular cushions (pair) size 1in, 1.5in & 2in - please ask for ITEM #'s