

foxdenton

Colour range

- Red - size 1
- Blue - size 2
- Green - size 3
- Beige - size 4
- Grey - size 5



Supportive and inclusive

Foxdenton is a solid classroom chair for children from the age of 3 years old who require minimal levels of postural support when seated. The appearance of the Foxdenton school chair has been purposely created to ensure your child feels more included in the classroom environment.

Effective learning

Playing, communicating, eating and carrying out desk-based activities are essential elements of school life. However, for some children these activities become so much harder when they are seated in a standard classroom chair. The Foxdenton provides support when they need it the most.

Adjustable support

Good seating posture starts with stabilising the pelvis. The Foxdenton has simple side supports that can be adjusted to keep your child's pelvis better aligned. A flip-up footrest also enables your child to sit in a more comfortable position.

Foxdenton technical data

all dimensions in mm unless otherwise stated

Code	Size	Seat height	Width	Seat depth	Weight limit (kg)
8011	1	300	120-240	255	20
8012	2	340	140-260	285	35
8013	3	380	180-310	320	45
8014	4	390	200-380	365	70
8015	5	430	200-380	365	80

foxdenton

accessories

The Foxdenton chair package comes complete with a school chair, Foxdenton seat with adjustable side supports, flip up footrest and a buckle pelvic strap.



Flip up footrest

CODE FR028 - FR030

Flips up to give clearance to facilitate standing transfer



Foot straps

CODE ST033 - ST034

Maintains foot position



Foxdenton seat

CODE FS001 - FS005

To facilitate the fitting of additional accessories



Adjustable side supports¹

CODE SP005

Provides hip and trunk support
Easy to adjust



Backrest¹

CODE BK137 - BK141 (40mm)

BK142 - BK146 (60mm)

Additional back support
Enables seat depth reduction



Pelvic strap

CODE ST035

Maintains hip position

Backrests technical data

all dimensions in mm unless otherwise stated

Size	1	2	3	4	5
40mm backrest	BK137	BK138	BK139	BK140	BK141
60mm backrest	BK142	BK143	BK144	BK145	BK146