

heathfield

The perfect start

Heathfield is the perfect starter chair for very young children with sizes up to young adults who require minimal levels of postural support when seated. Its unobtrusive design means that it fits well into nursery, school and home environments.

A natural choice

We believe in wood as a natural material for our chairs. Wood is extremely durable and it holds a contemporary beauty that makes it aesthetically pleasing while being environmentally friendly. Our timber comes from sustainable sources and offers thermal properties that make it warm to the touch.

Meeting your child's needs

The Heathfield and its many accessories work in partnership together to provide support and comfort for your child. We recognize that every child is different and the numerous configurations of the Heathfield will ensure that the chair is tailored to your child's needs.



The Heathfield chair comes complete with a buckle pelvic strap, adjustable armrests and a contoured seat.

Heathfield technical data

all dimensions in inches unless otherwise stated

Item #	Size	Seat to floor* height	Width	Seat depth	Back height	Max user weight (lbs)	Product weight (lbs)
6901	1	7.1 - 9.1	9.4	7.1 - 9.4	9.8	55.1	11
6902	2	7.9 - 11.8	11.0	8.7 - 12.2	11.0	99.2	13
6903	3	9.4 - 14.6	12.8	10.2 - 13.4	12.6	110.2	18
6904	4	11.4 - 15.4	14.4	12.6 - 15.7	14.4	132.3	35
6905	5	12.8 - 16.7	17.1	14.6 - 18.1	16.5	154.3	44

* Height of chair with contoured seat. The height with a flat seat is 1 inch less

heathfield accessories



Skis

ITEM # SK001 - SK005


Provides greater stability to avoid tippy-ness



Ski & full footrest

ITEM # SK023 - SK027

Provides foot support and offers greater stability



Mobile dolly base

ITEM # DL006 - DL010

Easy to move with swivel multi-directional casters



Reversible footboard

ITEM # FB006 - FB010

Foot support that prevents unwanted chair movement



Adjustable footrest

ITEM # FR005 - FR009

Easy to adjust footrest
The sides of the footrest help to keep feet within the chair



Minimum raise mobile base

ITEM # MB005 - MB007

Facilitates standing transfer
Easy to maneuver
Only available with sizes 3 - 5 chairs



4 Point chest harness

ITEM # HN001 - HN004



Pelvic harness

ITEM # HN039 - HN041



heathfield accessories



Pelvic cushions – double thickness (pair)¹

ITEM # CS050 - CS053, CS230-
CS234, CS077 - CS080

Decreases seat width
Available in various sizes



Pelvic cushions – single thickness (pair)²

ITEM # CS005 - CS008, CS014 -
CS017, CS023 - CS026

Provides hip and thigh support
Available in various sizes



Pommel⁵

ITEM # PL001 - PL002, PL013 -
PL014

Encourages hip abduction



Footplate & straps³

ITEM # FP001 - FP005

Ensures feet are firmly in position



Knee blocks^{3,5}

ITEM # KN001

Maintains knee position
Encourages hip abduction
Only available with sizes 3 - 5 chairs



Adjustable height push handle

ITEM # PH010 - PH014

Adjusts to achieve optimum push position. To be used in conjunction with a Mobile dolly base or Minimum raise mobile base



Tray⁴

ITEM # TR001 - TR005

Promotes trunk support
Additionally can be used for activities,
playing and feeding



Quick release tray⁶

ITEM # TR101 - TR104

The tray comes in white and only
available on chairs sizes 1-4

1 Available in sizes 1+1in, 1.5+1in & 2+1in

2 Available in size 1in, 1.5in & 2in

3 Only available in conjunction with a Mobile dolly base, Reversible footboard or Skis with half & full footrest

4 Heathfield chair sizes 1 & 2 must have a suitable stabiliser (footrest, mobile dolly base or skis) when attaching trays

5 The Pommel and Knee Blocks cannot be used together

6 Heathfield chair sizes 1 & 2 must have a suitable stabiliser e.g. footrest/dolly/skis in order to accommodate tray