

corner chair



Maximum adjustability

The Corner Chair offers high level of adjustability in the back, side angle, side pad, headrest andommel. This means that the chair can be positioned to suit your child's individual needs while also enabling the chair to easily grow with your child. Maximum user weight of 187lbs.

Floor level support

Our Corner Chairs are ideal for young children who need to develop sitting balance. During floor sitting some children are inclined to fall back or sideways if they overstretch or lose concentration. Our Corner Chairs are designed to provide support at the back and sides of children to improve their sitting balance. This encourages social interaction and allows them to engage in floor-based activities without the risk of falling over. Theommel at the front is designed to keep legs abducted to promote a stable posture.

Corner Chair technical data

all dimensions in inches

Item #	Size	Height range	Width	Seat depth
5501	1	16.1-20.4	15.7	11.8
5502	2	20.4-25.9	15.7	11.8
5503	3	21.7 - 27.2	17.7	13.7

hip spica

Support and comfort

Our range of Hip Spica chairs has been specifically designed for young children who are in plaster casts following hip surgery and require supportive seating.

Portable Hip Spica chair

The new design offers high levels of adjustability and creates a lightweight solution. The Portable Hip Spica features a new slimline design and can easily be dismantled so it's simple to transport. It is supplied with a 5-point harness and quick release tray as standard.



Portable Hip Spica chair technical data

all dimensions in inches

Item #	Size	Seat height	Footprint	Tray size	Max user weight lbs	Product weight lbs
5596	1	8.3 - 14.2	22.0 x 26.0	21.7 x 17.3	55	22