

juni

Colour range



plus natural wood (applies to tables and boosters only)



Posture at its foundation

Juni's seat and backrest both provide unique contouring. This promotes better posture, gives higher levels of comfort and helps your child feel more secure and stable in the Juni chair. It is our vibrant chair for children who require minimal levels of postural support when seated.

Highly adjustable

We understand that your child's needs change during each stage of growth and development. Seat height and depth as well as armrest height can all be easily adjusted. The backrest can be reclined by 10° (from 95° to 105°), while it also features a unique built-in lumbar support.

Dynamic learning

Whether in the classroom or at home, better posture whilst seated creates a more positive experience for your child; fine motor skills develop, confidence increases and the ability to engage and learn becomes easier.

Juni technical data

all dimensions in mm unless otherwise stated

Size	1	2	3	4
Height	160-220	200-280	240-340	300-400
Width	240	280	325	365
Depth	180-240	220-300	260-360	320-420
Back height	260	290	325	400
Back angle (°)	95/100/105	95/100/105	95/100/105	95/100/105
Weight limit (kg)	30	40	50	60
Colour options	Codes			
Red	16011	16012	16013	16014
Lime	16021	16022	16023	16024
Blue	16031	16032	16033	16034



Reclining backrest



Contoured seat



Height adjustable



Juni adjustable & foldable table

The Juni table offers two key benefits; the legs are adjustable so your child can work or play at the correct height and the table folds, which allows for easy storage. With a variety of colours (lime, blue, red or natural wood) this is the perfect companion to your Juni chair.

Juni table technical data

all dimensions in mm unless otherwise stated

Size	1	2
Height	310-450	400-600
Width/length	700	780
Depth	530	590
Colour options	Codes	
Red top	2411	2412
Lime top	2421	2422
Blue top	2431	2432
Natural top	2441	2442

juni accessories

A no fuss tray

The Quick release tray has been designed with the carer in mind. The tray is easy to attach and easy to remove. These two actions can all be accomplished by just using one hand, leaving your other hand free. It is also child proof, which means that your child cannot release the tray.

Our Quick release tray is light and easy to clean. Your child will be able to enjoy a variety of activities, such as playing and feeding. The tray also provides extra trunk support.



Quick release tray

Patent applied for
CODE TR101 - TR104

The tray is available in white

Juni Boosters technical data

all dimensions in mm unless otherwise stated

Booster Type	Hi Booster		Lo Booster	
Size of chair on booster ¹	1	2	1	2
Seat height	670-770		870-970	
Footprint open	520 x 520		520 x 640	
Folded size	520 x 340 x 700		520 x 420 x 930	
Colour options	Codes			
Red	JB112	JB122	JB111	JB121
Lime	JB212	JB222	JB211	JB221
Blue	JB312	JB322	JB311	JB321
Natural	JB412	JB422	JB411	JB421

¹ Please specify on your order size of the Juni chair that will be used on the Juni Booster

Colour range



plus natural wood

The new Juni Boosters

The Juni chair just got better with its new companions the Juni Boosters. Our Juni Boosters provide extra height for your child.

Hi Booster

Perfect for the breakfast bar or for feeding your child without crouching or kneeling.



Only available for Juni size 1 and 2

Lo Booster

Perfect for sitting your child at the dinner table or giving your child extra height.



Only available for Juni size 1 and 2



The Juni chair comes with contoured seat and backrest, built-in lumbar support, pelvic strap, height and depth adjustable seat, height adjustable armrests as well as angle adjustable backrest as standard.



Skis

CODE SK101 - SK104

Helps to stabilise and prevent chair movements



Mobile skis

CODE DL101 - DL104

Gives clearance to facilitate standing transfer
Easy to move



Ski & footboard¹

CODE SK111 - SK114

Provides foot support and offers greater stability



Sliding footrest²

CODE DL111 - DL114

Provides foot support
Slides backwards to give clearance to facilitate standing transfer



Height adjustable push handle

CODE PH101 - PH104

Adjusts to achieve optimum push position



Side support cushions³

Provides hip and thigh support
Available in 25mm and 50mm



Pommel³

Encourages hip abduction



4-point harness³

Provides trunk support