

# mobi changer



This image shows the product in a stainless steel option

## Room to room

The Mobi Changer is a fully mobile changing table that requires no installation or attachments. Its compact and slimline design means it can be easily maneuvered from room to room, offering an extremely cost effective option for multi-transfer or multi-user environments.

## Height adjustable

The electrically operated Mobi Changer provides an excellent travel range of 37.4in in height and can be lowered and elevated at the touch of a button. Safe and suitable working heights can be easily achieved which helps prevent the risk of back injury to the caregiver.

## Easy transfers

The Mobi Changer has been designed to ensure sufficient clearance for a wide number of hoists. Its mobility and unique stowaway side rails enable safe and easy handling. Whether you are transferring from a bed, wheelchair or toileting chair, the Mobi Changer is the perfect solution.

# mobi changer

## key features



This image shows the product in a powder coated option



- Mobile and easy to maneuver
- No installation required
- Cream powder coated frame for dry rooms
- Stainless steel for wet rooms
- Available in 2 sizes
- Length of table available 55.1in or 70.9in
- Weight capacity of 500.4lbs
- Travel range of 37.4in in height
- Easy to clean padded surfaces
- Optional angled back
- Unique stowaway side rails

Mild Steel Powder Coated	Stainless Steel	Size	User Weight	Table Dimensions		Travel Range	Safety Sides		Safety Sides Operational Space	Castor Centres	Nominal Clearances	Shipping Weights	
Item #	Item #		lbs	L in	W in	H in	Width X in	Height Y in	Z in	P in	Q in	G in	lbs
3121	3141	1	500.4	55.1	29.5	18.5 - 37.4	33.1	6.5	7.9	42.9	27.6	7.5	189.6
3123	3143	3	500.4	70.9	29.5	18.5 - 37.4	33.1	6.5	7.9	42.9	27.6	7.5	198.4

