

mobi heavy duty changer



Meeting bariatric needs

The Mobi Heavy Duty Changer has been designed to meet the needs of bariatric users. Its powerful 12K Newton motor, coupled with the extra wide mattress platform, makes it the ideal solution.

Maneuverability

The benefit of the Mobi range of products is that they do not require any wall mounting or floor attachments. They are easy to maneuver which means that the user can be transported from room to room with no trouble (door opening has to be 43.3in or more).

Height adjustable

The Mobi Heavy Duty Changer is electronically operated, which makes adjusting the height of the table simple, this is beneficial for both the caregiver and user.

mobi heavy duty changer

key features



- Mobile and easy to maneuver
- No installation required
- Electronically operated
- Cream powder coated frame for dry rooms only
- Length of table available in 78.7in
- Weight capacity of 1000lbs
- Travel range of 36.2in in height
- Easy to clean padded surfaces
- Unique stowaway side rails
- Wall mounted charger and 2 batteries come as standard

Mild Steel Powder Coated Item #	Size	User Weight lbs	Table Dimensions		Travel Range H in	Safety Sides		Safety Sides Operational Space Z in	Caster Centres		Nominal Clearances G in	Shipping Weights lbs
			L in	W in		Width X in	Height Y in		P in	Q in		
3161	1	1000	78.7	37.8	21.7 - 36.2	41.3	6.5	7.9	59.1	31.5	5.9	220.5

