

mobi shower trolley



Safe showering

For many children and adults, conventional methods of bathing and showering can be difficult and even hazardous. The Mobi Shower Trolley is a powered, height-adjustable shower trolley specifically designed to enable the user to be safely and easily washed or showered by their caregiver.

Reduces back injury

The Mobi Shower Trolley offers a range of height adjustments to enable the caregiver to shower and change the user at a safe and suitable working height. This reduces the risk of back injury for the caregiver.

Stowaway side rails

The Mobi Shower Trolley has been designed with two unique safety side rails that can be easily stowed away beneath the trolley. This allows maximum clearance for safe, easy transfer and use.

mobi shower trolley

key features



- Mobile and easy to maneuver
- No installation required
- Stainless steel frame
- Drainage channel and pipe
- A soft comfortable mattress and wedge shaped pillow
- Available in 1 size
- Length of table is 70.9in
- Weight capacity of 485lbs
- Travel range of 36.4in in height
- Easy to clean padded surfaces
- Unique stowaway side rails

Stainless Steel	Size	User Weight	Table Dimensions		Travel Range	Safety Sides		Safety Sides Operational Space	Caster Centres		Nominal Clearances	Shipping Weights
			L in	W in		Width X in	Height Y in		Z in	P in		
3151	1	485.0	70.9	29.5	18.3 - 36.4	33.1	7.7	7.9	42.9	27.6	7.5	189.6

