

whirl terzo

Colour range



Correct standing

Our new Whirl Terzo promotes independence but also encourages inclusion and helps children interact with their peers. The compact design makes the Whirl Terzo the perfect choice of stander for home and school. Suitable for children from 12 months old to young adults.

Innovative design

The aim with the Whirl Terzo was to create a superior upright stander at an affordable price. Reaching the optimum standing position for your child is key, the 10° tilt (which is easily controlled by a gas strut) and the high levels of adjustability will assist you in achieving this.



Improved manoeuvrability

The Whirl Terzo is an extremely compact stander and its design makes it highly manoeuvrable. Its small footprint means it fits up close to a table, once the tray is removed. This allows your child to participate in a range of activities.

Whirl terzo technical data

all dimensions in mm unless otherwise stated

Size	1	2
Code	8349	8350
User weight limit (kg)	60	80
Thoracic support (top) to foot plate height	*630-1020	*760-1270
Pelvic support (top) to foot plate height	*410-830	*560-1030
Knee support (top) to foot plate height	*130-420	*165-520
Foot plate to floor height (excluding sandals)	70	70
Max width (between thoracic supports)	300	460
Tray size (depth and width)	610 x 390	700 x 800
Angle (°)	0 to 10	0 to 10
Base frame	800 x 600	965 x 710
Shipping weight (including packaging) (kg)	38	40

*Minimum heights taken with footplate raised 50mm from floor

whirl

accessories

The Whirl and Whirl Terzo comes complete with anterior, pelvic and knee supports, cocoon sandals, activity tray and bowl with cover.

Modular parts such as the anterior supports, de-rotational pelvic supports and shoulder harness are interchangeable and work on both the Whirl and Whirl Terzo.



Knee supports¹

CODE KN006 - KN010

Height, width and depth adjustable to maintain knee position



De-rotational pelvic supports²

CODE SP003, SP004, SP022, SP023

Adjusts and maintains hip position to control obliquity and rotation into a midline position



Anterior supports²

CODE SP001, SP002, SP020, SP021

Adjustable to maintain trunk position



Shoulder harness²

CODE HN015 - HN018

Provides support from the rear to encourage symmetry and spinal extension



Tray grab rail full width

CODE GB001 - GB002

Helps to stabilise limbs and trunk



Suspension rail

CODE SR001 - SR002

Can be used to attach toys for stimulation

¹ Available in sizes 1 - 5

² Available in sizes 1 - 4