


zoomi

Benefits of Zoomi

- Economical seating solution
- Grows with your child
- Highly adjustable
- Tailored to the needs of your child
- Contemporary design
- Easy to clean
- Stable base with slimline footprint
- Can be used at home and in school

Our standard Zoomi

The standard Zoomi is ideally suited to the classroom and home. Both the seat and footrest can be easily adjusted to the appropriate height and depth, allowing the chair to accommodate your child's needs as they grow.



Standard features

- Frame
- Low backrest
- Flat laminated seat
- Laminated footboard
- Pelvic strap

Zoomi technical data

all dimensions in inches unless otherwise stated

Size	1	
Seat to floor height range	4.3 - 24.2	
Seat width	17.7	
Seat depth	21.7	
Back height	10.2	
Back angle (°)	95/100/105	
Frame height	31.9	
Max user weight (lbs)	187.4	
Product weight (lbs)	35	
Colour	Zoomi standard	Zoomi +
Blue	13031	13032

Perfectly customizable

With multiple options to cater for a variety of users, Zoomi can be uniquely configured to meet the needs of your child. We have put together two simple options:

- Standard Zoomi
- Zoomi+

zoomi+

One size fits all

Zoomi's innovative and unique seating system has been designed to grow with your child, covering an extensive age range, all the way up to young adults.

Zoomi+

Following collaboration with a number of therapists, we have developed Zoomi+. Complete with low back padded cover, armrest, quick release tray, 5-point harness, flat padded seat, slotted footrest and mobile skis, it offers the ideal solution for children who require a little extra postural support.

Standard features

- Frame
- Low back padded cover
- Armrest
- Quick release tray
- 5-point harness
- Flat padded seat
- Slotted footboard
- Mobile skis




Armrest

ITEM # AR019

Provides arm support

Used for the Quick release tray



Quick release tray

ITEM # TR104

Helps to stabilize trunk

Available in white only

Must be used with armrest